


 Trusted Clinical
Solutions



Biomaterials

Portfolio



 **ZimVie**

Trusted Clinical Solutions.



Table of Contents

THE POWER OF PUROS ALLOGRAFTS

Introduction	2-3
--------------------	-----

ALLOGRAFT BONE GRAFTS

Puros® Allografts	4
Cancellous Particulate Allograft	4
Cortical Particulate Allograft	5
Cortico-Cancellous Particulate Allograft	6
Puros DBM and Puros Ci Particulate Allograft	7
Puros Customized Blocks Bone Allograft	8
Puros Bone Block Allograft	9
Freeze-Dried Allograft Bone Grafts	10
RegenerOss® CC Cortico-Cancellous Particulate	10
RegenerOss® Particulate	11
DBM Putties	12
RegenaVate® Formable DBM	12
RegenerOss® Allograft Putty Plus	13
Xenogenic Bone Grafts	14
RegenerOss® Resorbable Xenograft Porcine Anorganic Bone Mineral	14
Endobon® Xenograft Bovine Granules	15
Combination Bone Grafts	16
RegenerOss® Bone Graft Plug	16
Synthetic Bone Grafts	17
IngeniOs® HA Synthetic Bone Particles	17
IngeniOs® β-TCP Bioactive Synthetic Bone Particles	18

SOFT TISSUE GRAFTS

Puros® Dermis Allograft Tissue Matrix	19
---	----

RESORBABLE BARRIER MEMBRANES

Puros® Pericardium Membrane	20
CopiOs® Pericardium Membrane	21
CopiOs® Extend Collagen Membrane	22
Socket Repair Membrane	23
OsseoGuard® and OsseoGuard Flex® Collagen Membranes	24
BioMend® and BioMend Extend™ Collagen Membranes	25

NON-RESORBABLE BARRIER MEMBRANES

OsseoGuard® PTFE Titanium Reinforced Membrane	26-27
OsseoGuard® Titanium Mesh High-Density PTFE Membrane	28
OsseoGuard® PTFE Textured and Non-Textured Membranes	28

WOUND DRESSINGS

Collagen Matrices: Plug, Tape, and Patch	29
--	----

SUTURES AND INSTRUMENTS

OsseoGuard® Non-Resorbable Sutures	30
Safescraper Twist Cortical Bone Collector	31

INTRUMENTS

Sinus Crestal Lift Instrument Kit	32
Sinus Lateral Lift Instrument Kit	32
Screw Fixation Instrument Kit	33

The Power of Puros® Allografts


Clinicians around the globe have counted on the Puros® family of allografts for hard tissue augmentation procedures for years.

The brand's renowned reputation is based on:*

- Predictable processing and configuration
- Clinical use in dentistry since 1999¹⁻³
- Collectively backed up by more than 400 scientific articles¹⁻⁵
- Supporting creation of healthy, vital bone⁶⁻⁹
- Predictable remodeling shown in human clinical studies¹⁰⁻¹⁵
- Ease of use and terminal sterilization¹⁶
- Quick hydration, five-year shelf life, and storage at room temperature¹⁶

More Studies Than Any Other Allograft⁵

Up to **23.7%** more vital bone formation with Puros Cancellous Particulate Allograft and Puros Cortical Particulate Allograft (1:1 ratio) – compared to freeze-dried allograft bone in sinus lift procedures.²⁰



Visual Comparison of Puros Cancellous Allograft to Natural Bone in SEM Image




Fig. 1

SEM images at 20,000x magnification of:
(a) Bone**
(b) Puros Cancellous Allograft


** Osteoclast-resorbed surface of human bone received unfixated, disinfected in 70% ethanol, air-dried, and rinsed in PBS.

The collagen fibrils are visible for Puros Cancellous Allograft following Tutoplast Processing and are similar to those seen in natural bone.¹⁹

*Claims referenced apply to Tutoplast processed grafts. ¹ Gambini A. et al. Chir Organi Mov (1999) 84:359-66. ² Rocci A. et al. Quintessence International, Edizione Italiana (1999) 15:373-380. ³ Semergidis T. et al. Int. J. Oral Maxillofac Surg (1999) 28:91. ⁴ Baldi D. et al. Implant Dent (2019) 28:472-477. ⁵ PubMed search (July 6th, 2020). ⁶ Tsao Y.P. et al. J Periodontol (2006) 77:416-25. ⁷ Leonetti J.A. et al. Implant Dent (2003) 12:217-226. ⁸ Keith J.D. et al. Int J Periodont Rest (2006) 26:321-327. ⁹ La Monaca G. et al. Case reports in dentistry (2019) 8, Article ID 6725351. ¹⁰ Froum S.J. et al. Int J Periodont Rest (2006) 26:543-51. ¹¹ Noubissi S.S. et al. J Oral Implantol (2005) 31:171-9. ¹² Block M.S. et al. J Am Dent Assoc (2002) 133:1631-1638. ¹³ Minichetti J.C. et al. J Oral Implantol (2004) 30:74-82. ¹⁴ Schmitt C.M. et al. Clin Oral Implants Res (2013) 24:576-85. ¹⁵ Soardi C.M. et al. Int J Oral Maxillofac Implants (2016) 31:352-8. ¹⁶ Puros Allograft IFU latest revision. ¹⁷ Data on File with RTI Surgical Inc. ¹⁸ Tadic D. et al. Biomaterials (2004) 25:987-94. ¹⁹ Ajami E, et al. J Oral Implantol (2023) 38: 169-180. ²⁰ Monje A. et al. Int J Oral Maxillofac Implants (2017) 32:121-127.

The Proprietary Tutoplast® Process

In 1969 the Tutoplast Tissue Sterilization Process was developed to sterilize and preserve tissue for implantation. More than 11 million implants have been sterilized through the Tutoplast Process with zero confirmed incidence of implant-associated infection.¹⁷



The Benefits of the Multi-Step Tutoplast Process

For allograft bone grafts, the process preserves the valuable bone mineral, collagen matrix, and tissue integrity¹⁸ while inactivating pathogens and gently removing unwanted materials, such as cells, antigens, and viruses¹⁷ – resulting in predictable, reliable, sterile, and safe tissue.¹⁷

*Images depict dermal processing

Puros® Cancellous Particulate Allograft

With a history of well-documented clinical results, Puros Cancellous is an easy-to-handle choice for predictable bone reconstruction and acts as an osteoconductive scaffold for new bone formation.¹⁻⁸

Clinical Evidence

- Up to 127% more vital bone formation compared to non-resorbable xenograft in sinus-lift procedures^{2,3,9}
- Newly formed vital bone after 3 to 5 months^{4,8,10} in extraction sockets
- 56% more graft-to-bone contact compared to non-resorbable xenograft³
- 9.7 mm vertical gain after 4 to 5 months when using Puros Allograft particulate with tenting screws¹¹
- Retains osteoconductive properties due to the preservation of the natural bone matrix collagen and mineral composition, trabecular pattern, and original porosity;^{1-6, 8, 12-14} enabling the ingrowth of vascular and cellular connective tissue⁴

Clinically successful in procedures for:

- Repair of periodontal bone and furcation defects^{1, 6, 15}
- Reconstruction of extraction sockets^{4, 7, 8, 10}
- Reconstruction of gaps around block grafts^{12, 13}
- Horizontal and vertical alveolar ridge augmentation¹⁶⁻¹⁹
- Sinus augmentation^{2, 9, 20, 21}



PUROS CANCELLOUS PARTICULATE ALLOGRAFT

Item Number	Description
68210	Puros Cancellous Particulate, 0.25 – 1 mm / 0.5 cc
68211	Puros Cancellous Particulate, 0.25 – 1 mm / 1 cc
68209	Puros Cancellous Particulate, 0.25 – 1 mm / 2 cc
68212	Puros Cancellous Particulate, 1 – 2 mm / 0.5 cc
68213	Puros Cancellous Particulate, 1 – 2 mm / 1 cc
68214	Puros Cancellous Particulate, 1 – 2 mm / 2 cc

Shelf-life: Five (5) years

¹ Tsao Y.P. et al. J Periodontol (2006) 77:416-25. ² Froum S.J. et al. Int J Periodontics Restorative Dent (2006) 26:543-51. ³ Noumbissi S.S. et al. J Oral Implantol (2005) 31:171-9. ⁴ Minichetti J.C. et al. J Oral Implantol (2004) 30:74-82. ⁵ Data on File with Rti Surgical Inc. ⁶ Dayi E. et al. J Int Med Res (2002) 30:168-73. ⁷ Baldi D. et al. Implant Dent (2019) 28:472-477. ⁸ Block M.S. et al. J Am Dent Assoc (2002) 133:1631-1638. ⁹ Schmitt C.M. et al. Clin Oral Implants Res (2013) 24:576-85. ¹⁰ Beck T.M. et al. J Periodontol (2010) 81:1765-72. ¹¹ Le B. et al. J Oral Maxillofac Surg (2010) 68:428-435. ¹² Keith J.D. et al. Int J Periodontics Restorative Dent (2006) 26:321-327. ¹³ Leonetti J.A. et al. Implant Dent. (2003) 12:217-226. ¹⁴ Tadic D. et al. Biomaterials (2004) 25:987-94. ¹⁵ Reddy B. et al. Journal of International Society of Preventive and Community Dentistry (2016) 6:248-253. ¹⁶ Block M.S. et al. J Oral Maxillofac Surg (2004) 62:67-72. ¹⁷ Le B. et al. Implant Dent (2008) 17:40-50. ¹⁸ Ronda M. et al. Clin Oral Implants Res (2014) 25:859-66. ¹⁹ La Monaca G. et al. Case reports in dentistry (2019) 8, Article ID 6725351. ²⁰ Soardi C.M. et al. Int J Periodontics Restorative Dent (2020) 40:757-764. ²¹ Monje A. et al. Int J Oral Maxillofac Implants (2017) 32:121-127.

Puros® Cortical Particulate Allograft

Puros Cortical can be used in space maintenance and volume enhancement procedures.^{1,2} It is slow-resorbing and maintains an open network for the proliferation of bone-forming cells.^{1,3}

Clinical Evidence

- Without sacrificing ridge contour, cortical particles remodel into a dense, lamellar structure as well as viable bone – with similar density to native bone⁴
- 2 mm in buccal bone thickness when used in a “sandwich” technique for the treatment of localized buccal dehiscence defects⁴
- 40% mineralized bone and 0.47% residual grafting materials after 4 months healing time in sinus lift procedures⁵
- Clinical and radiographic graft stability after 5 years follow up in sinus lift procedures⁶
- Reduced vertical and horizontal bone resorption when used in immediate implant placement extraction sites⁷

Clinically successful in procedures for:

- Sinus augmentation^{3,5,8,9}
- Alveolar ridge augmentation^{2,10,11}
- “Tenting” and “Sandwich” grafting techniques¹²⁻¹⁶
- Immediate implant post extraction sockets⁷

PUROS CORTICAL PARTICULATE ALLOGRAFT

Item Number	Description
68271	Puros Cortical Particulate, 0.5 – 1 mm / 0.5 cc
68272	Puros Cortical Particulate, 0.5 – 1 mm / 1 cc
68273	Puros Cortical Particulate, 0.5 – 1 mm / 2 cc
68274	Puros Cortical Particulate, 1 – 2 mm / 0.5 cc
68275	Puros Cortical Particulate, 1 – 2 mm / 1 cc
68276	Puros Cortical Particulate, 1 – 2 mm / 2 cc

Shelf-life: Five (5) years



¹Wang H.L. et al. *Implant Dent* (2006) 15:8-17. ²El Chaar E. et al. *Int J Periodontics Restorative Dent* (2019) 39:491-500. ³Berberi A. et al. *Journal of Maxillofacial and Oral Surgery* (2015) 14:624-629. ⁴Park S.H. et al. *Int J Periodont. Rest* (2006) 26:589-95. ⁵Berberi A. et al. *Implant Dent.* (2016) 25:353-60. ⁶Annibaldi S. et al. *Implant Dent* (2011) 20:445-54. ⁷Orti V. et al. *J Periodontal Implant Sci* (2016) 46:291-302. ⁸Soardi C.M. et al. *Int J Periodontics Restorative Dent* (2020) 40:757-764. ⁹Monje A. et al. *Int J Oral Maxillofac Implants* (2017) 32:121-127. ¹⁰Abed P.F. et al. *J Int Acad Periodontol* (2020) 22:11-20. ¹¹Wen S. et al. *Int J Periodontics Restorative Dent* (2018) 38:79. ¹²Leong D.J. et al. *Implant Dent* (2015) 24:4-12. ¹³Fu J.H. et al. *Clin Oral Implants Res* (2014) 25:458-67. ¹⁴Fu J.H. et al. *Clin Oral Implants Res* (2014) 26:1150-7. ¹⁵Fu J.-H. et al. *Clin Adv Periodontics* (2012) 2:172-177. ¹⁶Lee A. et al. *Implant Dent* (2009) 18:282-90.

Puros® Cortico-Cancellous Particulate Allograft

An anatomic-based mix of 70% cortical and 30% cancellous bone particulate. This mixture combines the clinical advantages of both Puros Cortical and Puros Cancellous Particulate Allograft materials.

Key Attributes

- Ideal for filling large and small volume bony voids
- Retains osteoconductive properties due to the preservation of the natural bone matrix collagen and mineral composition, trabecular pattern, and original porosity¹
- Pre-mixed formulation, no need to mix on site
- Easy handling – quick hydration, five-year shelf life, and room-temperature storage
- Bone from single donor⁴

Clinically successful in procedures for:

- Maxillary sinus floor augmentation¹
- Vertical ridge augmentation around dental implants²
- Immediate implant placement into extraction sockets³



PUROS CORTICO-CANCELLOUS PARTICULATE ALLOGRAFT

Item Number	Description
68800	Puros Cortico-Cancellous Particulate, 0.5 cc / 0.25 – 1 mm
68801	Puros Cortico-Cancellous Particulate, 1.0 cc / 0.25 – 1 mm
68802	Puros Cortico-Cancellous Particulate, 2.0 cc / 0.25 – 1 mm
68803	Puros Cortico-Cancellous Particulate, 0.5 cc / 1 – 2 mm
68804	Puros Cortico-Cancellous Particulate, 1.0 cc / 1 – 2 mm
68805	Puros Cortico-Cancellous Particulate, 2.0 cc / 1 – 2 mm

¹ Soardi et al. (2011). Atrophic maxillary floor augmentation by mineralized human bone allograft in sinuses of different size: an histologic and histomorphometric analysis. Clin Oral Implants Res 22, 560-566. ² Ronda et al. (2013). Expanded vs. dense polytetrafluoroethylene membranes in vertical ridge augmentation around dental implants: a prospective randomized controlled clinical trial. Clin Oral Implants Res 25, 859-866. ³ Sarnachiaro et al. (2016). Immediate Implant Placement into Extraction Sockets with Labial Plate Dehiscence Defects: A Clinical Case Series. Clin Implant Dent Relat Res 18, 821-829.

Puros® DBM and Puros® Ci Particulate Allograft

Puros Ci Particulate combines DBM Particulate with Cortico-Cancellous bone for an optimal blend of scaffold and osetoinductive potential.

Key Attributes

- Puros Ci is a combination of mineralized and demineralized allograft bone (DBM)
- Only DBM with verified osteoinductive (OI) potential is used in processing
- Mineralized component acts as a scaffold; mineralized particles include both cancellous and cortical bone chips
- Proprietary Tutoplast processing preserves the native collagen matrix of the mineralized particles
- Proprietary Cancellous SP® Sterilization Process preserves OI potential
- Single donor source for mineralized and demineralized particles
- Terminally sterilized to SAL 10⁻⁶
- Convenient handling: quick hydration and room temperature storage²

PUROS DBM PARTICULATE ALLOGRAFT

Item Number	Description
DBMPART025	Puros DBM Particulate Allograft, 0.25 cc
DBMPART050	Puros DBM Particulate Allograft, 0.5 cc
DBMPART100	Puros DBM Particulate Allograft, 1.0 cc

PUROS CI PARTICULATE ALLOGRAFT

Item Number	Description
69800	Puros Ci Particulate Allograft, 0.5 cc / 0.25 – 1 mm
69801	Puros Ci Particulate Allograft, 1.0 cc / 0.25 – 1 mm
69802	Puros Ci Particulate Allograft, 2.0 cc / 0.25 – 1 mm
69803	Puros Ci Particulate Allograft, 0.5 cc / 1 – 2 mm
69804	Puros Ci Particulate Allograft, 1.0 cc / 1 – 2 mm
69805	Puros Ci Particulate Allograft, 2.0 cc / 1 – 2 mm



*DBM is evaluated for bone formation utilizing an in-vivo athymic rat model. Findings from an animal model are not necessarily predictive of human clinical results.¹ Data on file with RTI Surgical, Inc. ² Refer to labeling for specific storage parameters.

Puros® Customized Blocks

Bone Allograft

Custom-made blocks of Tutoplast-processed cancellous bone are processed using CAD/CAM technology based on a CBCT/CT scan of the defect area. This makes the procedure more comfortable for your patient by reducing surgery time.¹

Clinical Evidence


- Customized block fits precisely and congruently to the defect²
- Large contact surface area improves ingrowth of blood vessels and revascularization³
- Additional manual adjustment of the defect and of the customized block is seldomly required, allowing for reduced surgery time and reduced morbidity⁴
- Clinical reports have shown stable bone levels up to 2 years follow-up after implant placement^{5,6}

Clinically successful in procedures for:

- Horizontal and vertical ridge reconstruction^{2,5,6}

PUROS ALLOGRAFT CUSTOMIZED BLOCK

Item Number	Description
68217	Puros Allograft Customized Block Standard, 27 x 15 x 15 mm (max)
68218	Puros Allograft Customized Block Large, 60 x 30 x 30 mm (max)



¹ Schlee M. et al. *Implant Dent* (2013) 22:212-8. ² Würzler K.K. et al. *Implantologie Journal* (2015) 5:30-36. ³ Mcallister B.S. et al. *J Periodontol* (2007) 78:377-96. ⁴ Parthasarathy J. *Ann Maxillofac Surg* (2014) 4:9-18. ⁵ Engler-Hamm D. *Implantologie* (2018) 26:231-242. ⁶ Blume O. et al. *J Esthet Restor Dent* (2018) 30:474-479.

Puros® Bone Block

Bone Allograft

By eliminating the need to excise an autogenous block graft, these prefabricated blocks may save time and shorten the patient's rehabilitation period.

Clinical Evidence

- Retains osteoconductive properties due to the preservation of the natural bone matrix collagen and mineral composition, trabecular pattern, and original porosity^{2,3}
- Implants can be placed 5 to 6 months after grafting^{2,4}
- Prospective studies showing comparable results to grafting with autogenous bone blocks^{1,5,6}
- Restores volume to severely resorbed ridges effectively as shown after 9 years follow up^{1,2,4,7}

Clinically successful in procedures for:

- Horizontal bone grafting^{1,2,8,9}
- Vertical bone grafting^{4,5}

PUROS BONE BLOCK ALLOGRAFT

Item Number	Description
68220	Puros Bone Block Allograft, 10 x 18 x 8 mm
68221	Puros Bone Block Allograft, 15 x 18 x 8 mm



¹Schlee M. et al. Head & Face Medicine (2014) 10:21. ²Keith J.D. et al. Int J Periodontics Restorative Dent (2006) 26:321-327. ³Tadic D. et al. Biomaterials (2004) 25:987-94. ⁴Leong D.J. et al. Implant Dent (2015) 24:4-12. ⁵Laino L. et al. Biomed Res Int (2014) 2014:982104. ⁶Motamedian S.R. et al. Ann Maxillofac Surg (2016) 6:78-90. ⁷Bauchet T. Implant (2020) 26:1-8. ⁸Jacotti M. et al. Implant Dent (2012) 21:444-8. ⁹Tresguerres F.G.F. et al. Clin Implant Dent Relat Res (2019) 21:1087-1098.

RegenerOss® CC Cortico-Cancellous Particulate

An anatomic-based mix of cortical and cancellous particulate that can be used to fill bony voids in a variety of dental applications. The graft has been sterilized using the validated Cancell SP® Sterilization Process, lyophilized, and terminally sterilized.

Key Attributes

- Graft has a natural donor-based mix of cortical and cancellous particles
- Sourced from a single donor
- The proprietary Cancell SP® Sterilization Process inactivates and/or removes bacteria, viruses, fungi, and spores without the use of antibiotics
- Designed to preserve biological integrity and natural collagen structure of bone
- Scientifically proven and clinically successful¹
- Low-dose gamma irradiation is applied terminally

Clinical applications:

- Reconstruction of extraction sockets
- Ridge augmentation and reconstruction
- Repair of periodontal and peri-implant defects
- Sinus floor elevation



REGENEROSS CC CORTICO-CANCELLOUS PARTICULATE

Item Number	Description
RMCCS050	Cortico-Cancellous Particulate, Small, 0.5 cc / 0.125 – 1.0 mm
RMCCS100	Cortico-Cancellous Particulate, Small, 1.0 cc / 0.125 – 1.0 mm
RMCCS200	Cortico-Cancellous Particulate, Small, 2.0 cc / 0.125 – 1.0 mm
RMCCLO50	Cortico-Cancellous Particulate, Large, 0.5 cc / 1.0 – 2.0 mm
RMCCLO100	Cortico-Cancellous Particulate, Large, 1.0 cc / 1.0 – 2.0 mm
RMCCLO200	Cortico-Cancellous Particulate, Large, 2.0 cc / 1.0 – 2.0 mm

¹ Data on File at ZimVie Dental.

RegenerOss®

Particulate Allograft

An allograft particulate that is aseptically-processed with a proprietary technique that removes unwanted cells, yet preserves valuable lipids. Available in a broad range of configurations for regenerating a variety of sites.

REGENEROSS CANCELLOUS PARTICULATE ALLOGRAFT

Item Number	Description
RMCA205	Cancellous, Mineralized, 0.5 cc / 200 – 300 µm
RMCA305	Cancellous, Mineralized, 0.5 cc / 300 – 500 µm
RMCA505	Cancellous, Mineralized, 0.5 cc / 500 – 800 µm
RMCA210	Cancellous, Mineralized, 1 cc / 200 – 300 µm
RMCA310	Cancellous, Mineralized, 1 cc / 300 – 500 µm
RMCA510	Cancellous, Mineralized, 1 cc / 500 – 800 µm
RMCA220	Cancellous, Mineralized, 2 cc / 200 – 300 µm
RMCA320	Cancellous, Mineralized, 2 cc / 300 – 500 µm
RMCA520	Cancellous, Mineralized, 2 cc / 500 – 800 µm

REGENEROSS CORTICAL PARTICULATE ALLOGRAFT

Item Number	Description
RMCO205	Cortical, Mineralized, 0.5 cc / 200 – 300 µm
RMCO305	Cortical, Mineralized, 0.5 cc / 300 – 500 µm
RMCO505	Cortical, Mineralized, 0.5 cc / 500 – 800 µm
RMCO210	Cortical, Mineralized, 1 cc / 200 – 300 µm
RMCO310	Cortical, Mineralized, 1 cc / 300 – 500 µm
RMCO510	Cortical, Mineralized, 1 cc / 500 – 800 µm
RMCO220	Cortical, Mineralized, 2 cc / 200 – 300 µm
RMCO320	Cortical, Mineralized, 2 cc / 300 – 500 µm
RMCO520	Cortical, Mineralized, 2 cc / 500 – 800 µm
RDCO205	Cortical, Partially Demineralized, 0.5 cc / 200 – 300 µm
RDCO305	Cortical, Partially Demineralized, 0.5 cc / 300 – 500 µm
RDCO505	Cortical, Partially Demineralized, 0.5 cc / 500 – 800 µm
RDCO210	Cortical, Partially Demineralized, 1 cc / 200 – 300 µm
RDCO310	Cortical, Partially Demineralized, 1 cc / 300 – 500 µm
RDCO510	Cortical, Partially Demineralized, 1 cc / 500 – 800 µm
RDCO220	Cortical, Partially Demineralized, 2 cc / 200 – 300 µm
RDCO320	Cortical, Partially Demineralized, 2 cc / 300 – 500 µm
RDCO520	Cortical, Partially Demineralized, 2 cc / 500 – 800 µm



RegenaVate®

Formable DBM

This allograft putty contains demineralized bone matrix (DBM) and mineralized cortical cancellous bone chips with a porcine gelatin carrier. The graft is available in two forms: Room Temperature (RT) and Frozen, to meet clinician preference.

Key Attributes

Induces bone formation and facilitates bone growth*

- The DBM is tested for osseointegrity* in a scientifically-proven in vivo rat assay
- Unique DBM provides handling flexibility
- Mineralized chips provide for osteoconductivity
- Clinician can control graft consistency: gel, paste, or putty

Clinically successful in procedures for:

- Filling extraction sockets
- Alveolar ridge augmentation
- Sinus floor elevation




ROOM TEMPERATURE

Item Number	Description
005301Z	RegenaVate Formable DBM, RT, 1 cc
005302Z	RegenaVate Formable DBM, RT, 2 cc

FROZEN

Item Number	Description
001504Z	RegenaVate Formable DBM, 1 cm x 1 cm x 0.5 cm, 0.5 cc
001505Z	RegenaVate Formable DBM, 1 cm x 2 cm x 0.5 cm, 1 cc
001510Z	RegenaVate Formable DBM, 1 cm x 4 cm x 0.5 m, 2 cc



*These implants were evaluated in a human clinical study and were shown to induce bone formation. Each lot is tested using the athymic nude rat assay to verify osteoinductivity potential.

RegenerOss® Allograft Putty Plus


Comprised of a mixture of DBM and mineralized allograft to facilitate the regeneration of new bone. The plant-based carrier leaves no residual soy proteins.

Key Attributes

- Encourages bone growth potential by incorporating a high bone-to-carrier ratio without sacrificing its handling characteristics¹
- Contains 48% bone graft material by weight (28% cortical DBM and 20% cancellous mineralized bone chips)
- Moldable, non-toxic, lecithin carrier that is highly resistant to irrigation
- Osseointegrity of every lot is validated by a cell proliferation assay
- Ergonomic design features a smaller diameter syringe with a curved tip to treat hard-to-reach defects


Indications

- Pre-implant defects
- Extraction sockets
- Localized ridge augmentation
- Sinus augmentation



REGENEROSS ALLOGRAFT PUTTY PLUS

Item Number	Description
ROAPM05	RegenerOss Allograft Putty Plus, 0.5 cc Syringe
ROAPM10	RegenerOss Allograft Putty Plus, 1 cc Syringe
ROAPM20	RegenerOss Allograft Putty Plus, 2 cc (Four 0.5 cc syringes) Value Pack



¹ Urist, MR. Bone: Formation By Autoinduction. Science, 1965, 150: pp 893.

RegenerOss® Resorbable Xenograft

Porcine Anorganic Bone Mineral

RegenerOss Resorbable Xenograft has up to 95% porosity¹ enabling excellent osteoconductivity and adequate space for new bone formation.

Clinical Evidence

- Osteoconductive surface and interconnecting macro and microscopic porous structure that support the formation and ingrowth of new bone at the implantation site^{2,3}
- Clinical results showing new bone formation, both around and within the particles⁴
- Porcine-derived carbonate apatite shows superior osteoconductive potential than hydroxyapatite^{5,6}
- Resorption and remodelling profiles are closer to human bone than those of synthetic bone graft substitutes⁶

Clinically successful in procedures for:

- Augmentation around implants³
- Alveolar ridge augmentation/reconstruction^{3,7}
- Sinus lifts³
- Extraction sockets^{4,8-10}
- Periodontal defects³



REGENEROSS RESORBABLE XENOGRAFT

Item Number	Particle Size	Description
ROXR05	Small Particles	RegenerOss Resorbable Xenograft, 0.25 – 1 mm / 0.5 cc
ROXR10	Small Particles	RegenerOss Resorbable Xenograft, 0.25 – 1 mm / 1.0 cc
ROXR20	Small Particles	RegenerOss Resorbable Xenograft, 0.25 – 1 mm / 2.0 cc
ROXR40	Small Particles	RegenerOss Resorbable Xenograft, 0.25 – 1 mm / 4.0 cc
ROXRLG10	Large Particles	RegenerOss Resorbable Xenograft, 1 – 2 mm / 1.0 cc
ROXRLG20	Large Particles	RegenerOss Resorbable Xenograft, 1 – 2 mm / 2.0 cc
ROXRS025	Small Particles	RegenerOss Resorbable Xenograft, Syringe, 0.5 – 1 mm / 0.25 cc
ROXRS05	Small Particles	RegenerOss Resorbable Xenograft, Syringe, 0.5 – 1 mm / 0.5 cc

Shelf-Life Small & Large Particles: Three (3) years

Shelf-Life Syringe: Two (2) years

¹ Data on File with Collagen Matrix Inc. ² Klenke F.M. et al. J Biomed Mater Res A (2008) 85A:777-786. ³ RegenerOss Xenograft IFU latest revision. ⁴ Guarnieri R. et al. Regen Biomater (2017) 4:125-128. ⁵ Spence G. et al. J Biomed Mater Res A (2009) 90A:217-224. ⁶ Ellies L.G. et al. J Biomed Mater Res (1988) 22:137-48. ⁷ Cucchi A. et al. J. Oral Implantol. (2019) 45:59-64. ⁸ Guarnieri R. et al. J Oral Maxillofac Res (2019) 10:e3. ⁹ Guarnieri R. et al. J Oral Maxillofac Res (2017) 8:e5. ¹⁰ Lai V.J. et al. J Periodontol (2020) 91:361-368.

Endobon Xenograft

Bovine Granules

An essentially non-resorbable material that is ideally suited for regeneration of defects when effective space maintenance is required.¹

Clinical Evidence

- Fully deproteinized bovine-derived hydroxyapatite²
- Non-resorbable for predictable volume stability and maintenance³
- Using this in a buccal onlay tunnel technique showed 2-year ridge width after restoration was 9.8 +/-1.2 mm (range, 8.0 – 11.2 mm).¹
- Xenograft particles will be surrounded by newly formed vital bone⁴

Clinically successful in procedures for:

- Alveolar ridge augmentation, including aesthetic contouring defects^{1,5,6}
- Extraction socket grafting⁷
- Sinus elevation^{4,8}



ENDO BON XENO GRAFT GRANULES

Item Number	Particle Size	Description
ROX05	Small Particles	Endobon Xenograft Granules, 0.5 – 1 mm, 0.5 ml
ROX10	Small Particles	Endobon Xenograft Granules, 0.5 – 1 mm, 1 ml
ROX20	Small Particles	Endobon Xenograft Granules, 0.5 – 1 mm, 2 ml
ROXLG20	Large Particles	Endobon Xenograft Granules, 1 – 2 mm, 2 ml
ROXLG50	Large Particles	Endobon Xenograft Granules, 1 – 2 mm, 5 ml (5 units @ 1 ml each)
ROXLG80	Large Particles	Endobon Xenograft Granules, 1 – 2 mm, 8 ml (8 units @ 1 ml each)

Shelf-life: 18 months

¹ Block M.S. et al. J Oral Maxillofac Surg (2013) 71:1513-1519. ² Tadic D. et al. Biomaterials (2004) 25:987-94. ³ Block M.S. et al. J Oral Maxillofac Surg (2012) 70:1321-1330. ⁴ Nevins M. et al. Int J Periodontics Restorative Dent (2011) 31:227-35. ⁵ Barone A. et al. Int J Periodontics Restorative Dent (2013) 33:795-802. ⁶ Castillo R.a.D. Inside Dent (2011) 7:94-96. ⁷ Fischer K.R. et al. Int J Periodontics Restorative Dent (2018) 38:549-556. ⁸ Testori T. et al. Int J Periodontics Restorative Dent (2012) 32:295-301.

RegenerOss®

Bone Graft Plug

An easy-to-use grafting solution for filling extraction sockets and periodontal defects. The Plug features a combination of 80% graft particulate and 20% Type I bovine collagen that adapts to the shape of the defect once hydrated.

Key Attributes

- Plug-shaped form mimics the shape of extraction sockets for easy placement
- Combines 80% carbonate apatite granules with 20% Type I collagen from bovine Achilles tendon
- Mineral component of the graft has macro- and microstructures similar to human bone
- Collagen holds graft particles in place within the defect site
- Plug can be placed dry and takes the shape of the defect once hydrated

Clinical Applications

- Reconstruction of extraction sockets
- Ridge augmentation
- Repair of periodontal and peri-implant defects



REGENEROSS BONE GRAFT PLUG

Item Number	Description
RGP0625	RegenerOss Bone Graft Plug, 6 x 25 mm, Box of 5
RGP1020	RegenerOss Bone Graft Plug, 10 x 20 mm, Box of 5

Shelf-life: Three (3) years



IngeniOs HA

Synthetic Bone Particles

Long-lasting bone particles made of pure-phase hydroxyapatite (HA), a composition similar to HA found in naturally-occurring bone.¹

Clinical Evidence

- Significantly higher cell attachment was seen with IngeniOs HA compared to Geistlich Bio-Oss at all time points in an in-vitro study²
- Long-lasting osteoconductive support with negligible resorption to provide long-term graft stability and maintenance of volume
- Up to 80% interconnected porosity allowing for vascularized bone formation, osseointegration, and the natural remodeling process to occur within the graft framework^{3,4}
- Radiopacity of material making it easy to identify on an x-ray⁴
- Can be used as a bone graft extender to provide radiopacity or long-term volume preservation

Clinically successful in procedures for:

- Alveolar ridge augmentation/reconstruction^{1,4}
- Sinus lifts^{1,4}
- Defects after removal of bone cysts^{1,4}
- Extraction sockets⁴



INGENIOS HA SYNTHETIC BONE PARTICLES

Item Number	Description
0-802501	IngeniOs HA Synthetic Bone Particles, 0.25 – 1 mm / 0.25 cc
0-800501	IngeniOs HA Synthetic Bone Particles, 0.5 – 1 mm / 0.5 cc
0-801001	IngeniOs HA Synthetic Bone Particles, 0.5 – 1 mm / 1 cc
0-802001	IngeniOs HA Synthetic Bone Particles, 0.5 – 1 mm / 2 cc
0-900501	IngeniOs HA Synthetic Bone Particles, 1 – 2 mm / 0.5 cc
0-901001	IngeniOs HA Synthetic Bone Particles, 1 – 2 mm / 1 cc
0-902001	IngeniOs HA Synthetic Bone Particles, 1 – 2 mm / 2 cc

Shelf-life: Five (5) years

¹ Holweg A. et al. EDI Journal (2012) 3:64-73. ² Bernhardt A. et al. Clin Oral Implants Res (2011) 22:651-7. ³ Data on File with Curasan Ag. ⁴ IngeniOs HA Synthetic Bone Particles IFU latest revision.

IngeniOs β -TCP

Bioactive Synthetic Bone Particles

Resorbable bone particles made of pure-phase beta tricalcium phosphate (β -TCP) that is silicated, providing the potential for increased bioactivity.¹⁻³

Clinical Evidence

- Resorption and remodeling profiles of carbonate apatite mimic natural bone mineral¹
- Higher osteoconductive properties and earlier bioresorption, compared to HA samples^{2,3,4}
- Fully resorbable material designed to resorb in balance with replacement by naturally-regenerating mineralized bone³
- Up to 75% interconnected porosity to enable ingrowth of healthy bone tissue^{1,3}
- Radiopacity of material making it easy to identify on an x-ray³
- Can be used as a bone graft extender to extend volume or add radiopacity

Clinically successful in procedures for:

- Alveolar ridge augmentation/reconstruction³
- Filling of defects after root resection, apicectomy, and cystectomy⁶
- Sinus lifts³
- Extraction sockets³
- Periodontal defects³



INGENIOS β -TCP BIOACTIVE SYNTHETIC BONE PARTICLES

Item Number	Description
0-602501	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 0.25 – 1 mm, 0.25 cc
0-600501	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 0.5 – 1 mm, 0.5 cc
0-601001	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 0.5 – 1 mm, 1 cc
0-602001	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 0.5 – 1 mm, 2 cc
0-700501	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 1 – 2 mm, 0.5 cc
0-701001	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 1 – 2 mm, 1 cc
0-702001	IngeniOs β -TCP Bioactive Synthetic Bone Particles, 1 – 2 mm, 2 cc

Shelf-life: Five (5) years

¹ Data on File with Curasan Ag. ² Pietak A.M. et al. Biomaterials (2007) 28:4023-32. ³ IngeniOs β -TCP Bioactive Synthetic Bone Particles IFU latest revision.

Puros® Dermis

Allograft Tissue Matrix

A high-quality, natural, biocompatible dermal matrix used in horizontal and vertical soft-tissue augmentation.¹⁻³

Clinical Evidence

- After 5 years follow-up, no statistical significant differences in tissue thickening and gain of clinical attachment level compared to autogenous connective tissue graft when used to treat multiple gingival recessions¹
- Superior tissue characteristics due to solvent dehydration processing compared to freeze-dried grafts⁴
- Not cross-linked compared to a xenogeneic soft-tissue graft⁵
- 100% free of antibiotics: Puros Dermis tissue matrix is not treated with antibiotics like a certain freeze dried human dermis graft⁶
- Rehydration in a single bath reduces preparation time⁷

Clinically successful in procedures for:


- Horizontal and vertical soft-tissue augmentation^{1-3,8}
- Periodontal and peri-implant soft tissue management⁹⁻¹³

PUROS DERMIS ALLOGRAFT TISSUE MATRIX - THIN

Item Number	Description - Thin
68794	Puros Dermis Tissue Matrix, 10 x 10 mm, 0.3 – 0.8 mm
68795	Puros Dermis Tissue Matrix, 10 x 20 mm, 0.3 – 0.8 mm
68796	Puros Dermis Tissue Matrix, 10 x 40 mm, 0.3 – 0.8 mm
68797	Puros Dermis Tissue Matrix, 20 x 40 mm, 0.3 – 0.8 mm

PUROS DERMIS ALLOGRAFT TISSUE MATRIX - THICK

Item Number	Description - Thick
68793	Puros Dermis Tissue Matrix, 10 x 10 mm, 0.8 – 1.8 mm
68790	Puros Dermis Tissue Matrix, 10 x 20 mm, 0.8 – 1.8 mm
68791	Puros Dermis Tissue Matrix, 10 x 40 mm, 0.8 – 1.8 mm
68792	Puros Dermis Tissue Matrix, 20 x 40 mm, 0.8 – 1.8 mm



¹ Kroiss S. et al. Quintessence Int. (2019) 50:278-285. ² Petrungaro P. Inside Dent (2007) 3:2-4. ³ Petrungaro P.S. Inside Dent (2010) 2-9. ⁴ Hinton R. et al. Am J Sports Med (1992) 20:607-12. ⁵ Geistlich Fibro-Gide® IFU 08/2017. ⁶ Alloderm IFU 11/2017. ⁷ Puros Dermis Allograft Tissue Matrix IFU 06/2017. ⁸ Abou-Arraj R.V. et al. Int J Periodontics Restorative Dent (2017) 37:571-579. ⁹ Aroni M.a.T. et al. Rev Odontol UNESP (2016) 45:78-84. ¹⁰ Wang H.L. et al. J Periodontol (2014) 85:1693-701. ¹¹ Alasmari D.S. J Am Sci (2014) 10:97-99. ¹² Farina V. et al. Int J Oral Maxillofac Implants (2015) 30:909-17. ¹³ Puisys A. et al. Clin Oral Implants Res (2015) 26:123-9.

Puros® Pericardium Allograft Membrane

Allograft membrane provides a long-lasting barrier when an optimum balance of strength and handling for graft containment are necessary.^{1, 2, 4, 6}

Clinical Evidence

- Functions as a barrier during the critical part of wound healing and helps stabilize and maintain bone growth material in the defect space^{1, 2, 4, 6}
- Retains the natural collagen matrix and mechanical properties of native pericardium due to the proprietary Tutoplast Process
- Exhibits multi-directional strength
- Rehydrates quickly
- Three convenient sizes can be cut to shape for specific procedures
- Drapeable membrane adapts to the defect or grafted site

Clinically successful in procedures for:

- Guided bone regeneration procedures^{3, 4}
- General surgery applications⁷



PUROS PERICARDIUM MEMBRANE

Item Number	Description
68770	Puros Pericardium Allograft Membrane, 15 x 20 mm
68771	Puros Pericardium Allograft Membrane, 20 x 30 mm
68772	Puros Pericardium Allograft Membrane, 30 x 40 mm

Shelf-life: Five (5) years

¹ Sohn DS, Shin HI, Ahn MR, Lee JS. Piezoelectric vertical bone augmentation using the sandwich technique in an atrophic mandible and histomorphometric analysis of mineral allografts: a case report series. *Int J Periodontics Restorative Dent.* 2010;30(4):383-391. ² Taskonak B, Ozkan Y. An alveolar bone augmentation technique to improve esthetics in anterior ceramic FPDs: a clinical report. *J Prosthodont.* 2006;15(1):32-36. ³ Petrungero PS, Amar S. Localized ridge augmentation with allogenic block grafts prior to implant placement: case reports and histologic evaluations. *Implant Dent.* 2005;14(2):139-148. ⁴ Rocci A, Martignoni M. Local enlargement of the alveolar ridge using a mineralized allogenic cortical- cancellous block graft: a clinical case study. *Quintessence Int.* 1999;11(12):373-380. (Italian Edition). ⁵ Paolantonio M. Combined periodontal regenerative technique in human intrabody defects by collagen membranes and anorganic bovine bone. A controlled clinical study. *J Periodontol.* 2002 Feb;73(2):158-166. ⁶ Shin HI, Sohn DS. A method of sealing perforated sinus membrane and histologic finding of bone substitutes: a case report. *Implant Dent.* 2005;14(4): 328-335. ⁷ Keith JD, Salama MA. Ridge preservation and augmentation using regenerative materials to enhance implant predictability and esthetics. *Compend Contin Educ Dent.* 2007 Nov;28(11):614-621; quiz 622-624.

CopiOs® Pericardium Xenograft Membrane

A long-lasting, conformable barrier – strong enough to meet most clinical needs and supple enough to adapt to challenging graft contours.¹⁻⁴

Key Attributes

- Processed from bovine pericardium⁵
- Barrier time 8–24 weeks: for longer graft protection and stabilization^{1,6,7}
- Not side specific for convenient handling⁸
- Retains the structure and composition of natural pericardial tissue due to the proprietary Tutoplast Process^{9,10}
- High tensile strength and suture pull-out force may be useful in guided bone regeneration techniques⁷
- Clinically demonstrated performance in guided bone regeneration procedures where ease of manipulation and adaptability to surface contours is essential¹¹⁻¹⁴
- Shown to provide a stable, long-lasting barrier during healing and integration of bone graft materials, and staged or immediately placed implants^{12,15-17}
- Significantly thicker buccal bone plate when using CopiOs Pericardium membranes to cover bone graft during implant placement^{12,18}

Clinically successful in procedures for:

- Guided tissue regeneration (GTR) in periodontology^{5,19}
- Covering and protecting bone graft material, e.g. in guided bone regeneration procedures (GBR)^{5,11,12}

COPIOS PERICARDIUM MEMBRANE

Item Number	Description
77776	CopiOs Pericardium Membrane, 15 x 20 mm
77777	CopiOs Pericardium Membrane, 20 x 30 mm
77778	CopiOs Pericardium Membrane, 30 x 40 mm



¹ Rothamel D. et al. Clin Oral Implants Res (2005) 16:369-78. ² Data on file with RTI Biologics Inc, USA. ³ Leong D.J. et al. Implant Dent (2015) 24:4-12. ⁴ Berberi A. et al. J Maxillofac Oral Surg (2015) 14:263-70. ⁵ CopiOs Pericardium Membrane IFU latest revision. ⁶ Siar C.H. et al. Clin Oral Implants Res (2011) 22:113-20. ⁷ Gasser A. et al., Mechanical stability of collagen membranes: an in vitro study, in AADR/CADR Meeting, 2016: Los Angeles. ⁸ Data on File with Zimmer Biomet Dental. ⁹ Marashdeh M.Q.M., Characterization and Development of Optimization Strategy for the Processing of Allogenic and Xenogenic Bone and Pericardium. 2007, Thesis, University of Erlangen-Nürnberg. ¹⁰ Kasaj A. et al. Head Face Med (2008) 4:22. ¹¹ El Chaar E. et al. J Oral Implantol (2017) 43:114-124. ¹² Fu J.H. et al. Clin Oral Implants Res (2014) 25:458-67. ¹³ Soardi C.M. et al. Clin Adv Periodontics (2013) 4:1-7. ¹⁴ Fu J.-H. et al. Clin Adv Periodontics (2012) 2:172-177. ¹⁵ Sterio T.W. et al. Int J Periodontics Restorative Dent (2013) 33:499-507. ¹⁶ Le B. et al. J Oral Maxillofac Surg (2016) 74:1552-61. ¹⁷ Laino L. et al. Biomed Res Int (2014) 2014:982104. ¹⁸ Garaicoa C. et al. Clin Implant Dent Relat Res (2015) 17:717-23. ¹⁹ Schlee M. et al. Head Face Med (2012) 8:6.

CopiOs® Extend Collagen Membrane

CopiOs Extend Membrane is a long-lasting, resorbable collagen membrane designed to allow implant placement while providing ample time for regeneration.

Clinical Evidence

- Made of highly purified porcine dermis¹
- Barrier time 6 – 9 months¹
- Not side specific for convenient handling¹
- Cell-occlusive – allows nutrients to permeate while occluding epithelial cells²
- Convenient handling – conformable and easy to reposition in the defect
- Performs when primary closure has not been achieved³

Clinically successful in procedures for:

- Augmentation around implants placed in immediate and delayed extraction sockets¹
- Localized ridge augmentation for later implantation¹
- Alveolar ridge reconstruction for prosthetic treatment¹
- Filling of bone defects¹
- Guided bone regeneration in dehiscence defects¹
- Guided tissue regeneration procedures in periodontal defects¹



COPIOS EXTEND MEMBRANE

Item Number	Description
0190Z	CopiOs Extend Membrane, 15 x 20 mm
0191Z	CopiOs Extend Membrane, 20 x 30 mm
0192Z	CopiOs Extend Membrane, 30 x 40 mm

¹CopiOs Extend Membrane IFU latest revision. ²Data on File with Collagen Matrix Inc. ³Data on File with Zimmer Biomet Dental.

Socket Repair Membrane


Designed to assist wound healing in alveolar facial plate repair following atraumatic, flapless single-root tooth extraction.

Clinical Evidence

- Made of bovine achilles tendon¹
- Barrier time 26-38 weeks¹ (accelerated resorption will occur if exposed)
- Flapless approach preserves marginal soft-tissue contours² and does not compromise buccal bone tissue. Maintaining this tissue and the vascular supply to the area is important for achieving highly esthetic results³

Clinically successful in procedures for:

- 3-wall extraction sockets³⁻⁵



SOCKET REPAIR MEMBRANE

Item Number	Description
0154	Zimmer Socket Repair Membrane, 10 x 20 mm

¹ Zimmer Socket Repair Membrane IFU latest revision ² Danesh-Meyer M. Australasian Dental Practice (2008) 150-158. ³ Elian N. et al. Pract Proced Aesthet Dent (2007) 19:99-104. ⁴ Eskow A.J. et al. J Periodontol (2014) 85:514-24. ⁵ Hoang T.N. et al. J Periodontol (2012) 83:174-81.

OsseoGuard® and OsseoGuard Flex®

Collagen Membranes

Two levels of drapability for ease of use in various clinical procedures.

Clinical Evidence

- Made of bovine Achilles tendon (OsseoGuard)¹ and highly purified bovine dermis (OsseoGuard Flex)²
- Barrier time 6 – 9 months¹⁻³
- Not side specific for convenient handling⁴
- Can be trimmed, placed dry, or hydrated and finally sutured in place^{1, 2}
- Performs when primary closure has not been achieved (OsseoGuard Flex)⁴
- Space maintaining (OsseoGuard)⁵

OsseoGuard clinically successful in procedures for:

- Periodontal and/or dental surgery procedures¹
- In the area of periodontal defects, dental implant, bone defect or ridge reconstruction^{1, 6-9}

OsseoGuard Flex clinically successful in procedures for:

- Augmentation around implants placed in immediate extraction sockets, delayed extraction sockets^{2, 10-12}
- Localized ridge augmentation for later implantation^{2, 13}
- Alveolar ridge reconstruction for prosthetic treatment²
- Filling of bone defects²
- Guided bone regeneration in dehiscence defects²

OSSEOGUARD MEMBRANE


Item Number	Description
OG1520	OsseoGuard Resorbable Collagen Membrane, 15 x 20 mm
OG2030	OsseoGuard Resorbable Collagen Membrane, 20 x 30 mm
OG3040	OsseoGuard Resorbable Collagen Membrane, 30 x 40 mm

OSSEOGUARD FLEX MEMBRANE

Item Number	Description
OGF1520	OsseoGuard Flex Resorbable Collagen Membrane, 15 x 20 mm
OGF2030	OsseoGuard Flex Resorbable Collagen Membrane, 20 x 30 mm
OGF3040	OsseoGuard Flex Resorbable Collagen Membrane, 30 x 40 mm



OsseoGuard



OsseoGuard Flex

¹ OsseoGuard Membrane IFU latest revision. ² OsseoGuard Flex Membrane IFU latest revision. ³ Data on File with Collagen Matrix Inc. ⁴ Data on File with ZimVie. ⁵ Block M.S. et al. J Oral Maxillofac Surg (2013) 71:1513-1519. ⁶ Fischer K.R. et al. Int J Periodontics Restorative Dent (2018) 38:549-556. ⁷ Tan-Chu J.H. et al. Int J Periodontics Restorative Dent (2014) 34:399-403. ⁸ Block M.S. et al. J. Oral Maxillofac. Surg. (2012) 70:1321-1330. ⁹ Nevins M. et al. Int J Periodontics Restorative Dent (2011) 31:227-35. ¹⁰ Chasioti E. et al. Case reports in dentistry (2015) Article ID 439706:8pages. ¹¹ Castillo R.a.D. Inside Dent (2011) 7:94-96. ¹² Felice P. et al. Eur J Oral Implantol (2015) 8:375-84. ¹³ Chasioti E. et al. Quintessence Int (2013) 44:763-71.

BioMend® and BioMend Extend™ Collagen Membranes

Resorbable collagen membranes made of bovine Achilles tendon that are rigid enough to create and maintain space.¹

Clinical Evidence

- Two different options of barrier time: 8 weeks max. (BioMend), 18 weeks max. (BioMend Extend)²
- Not side specific for convenient handling³
- Cell-occlusive – serves as barrier to prevent epithelial cell migration and allows passage of essential nutrients²
- Up to 54% more horizontal bone gain when using BioMend Extend membranes to cover bone graft during implant placement⁴

*Compared to a Porcine Membrane**

- Significantly higher tensile strength in wet and dry state may be useful for guided bone regeneration techniques⁹
- 34% more new bone fill and 28% more bone to-implant contact when using BioMend Extend Membranes for treatment of implant dehiscence defects¹

Clinically successful in procedures for:

- Guided tissue regeneration procedures in periodontal defects²
- Periodontal surgery^{2,5,6}
- Use in dental surgery procedures as a material for placement in the area of an implant, bone defect, or ridge construction^{2,7}
- Sinus lift procedures⁸

BIOMEND MEMBRANE

Item Number	Description
0103Z	BioMend Resorbable Collagen Membrane, 15 x 20 mm
0105Z	BioMend Resorbable Collagen Membrane, 20 x 30 mm
0107Z	BioMend Resorbable Collagen Membrane, 30 x 40 mm

BIOMEND EXTEND MEMBRANE

Item Number	Description
0140Z	BioMend Extend Resorbable Collagen Membrane, 15 x 20 mm
0141Z	BioMend Extend Resorbable Collagen Membrane, 20 x 30 mm
0142Z	BioMend Extend Resorbable Collagen Membrane, 30 x 40 mm



¹Bio-Gide Membrane, Edward Geistlich Sohne AG. ²Oh T.J. et al. Clin Oral Implants Res (2003) 14:80-90. ³BioMend and BioMend Extend Absorbable Collagen Membrane IFU latest revision. ⁴Data on File with Collagen Matrix Inc. ⁵Park S.H. et al. Clin Oral Implants Res (2008) 19:32-41. ⁶Wang H.L. et al. J Periodontol (1994) 65:1029-36. ⁷Wang H.-L. et al. Periodontol 2000 (2012) 59:140-157. ⁸Saravanan P. et al. J Oral Implantol (2013) 39:455-62. ⁹Ranaan J. et al. Clin Oral Implants Res (2018). ⁹Coic M. et al. Rev Stomatol Chir Maxillofac Chir Orale (2010) 111:286-290.

OsseoGuard® PTFE Titanium Reinforced Membrane

OsseoGuard PTFE Titanium Reinforced High-Density PTFE Membranes are designed for periodontal applications, large defects, and defects missing adequate bony architecture.

The traditional frame design, incorporating delicate and strategically-placed titanium “struts”, has more than 25 years of clinical history and successful use in GBR. This innovative, hybrid design consists of a thin layer of expanded PTFE (e-PTFE) laminated to a textured d-PTFE membrane. In between these two layers lies a titanium framework. The titanium framework is a grade of titanium that has little to no memory. Once formed, the titanium-reinforced membrane will remain in that shape until mechanically altered.

Key Attributes¹


- Less Titanium Bulk
- Grade 1 titanium, lightweight framework
- Textured surface, d-PTFE backing
- Can be molded and shaped for tenting and space maintenance
- May be easily cut with scissors to custom-fit various defects

Benefits¹

- Less is more - less titanium bulk allows for greater versatility in shaping and placement
- Easy to form in three dimensions and retains no memory, allowing for passive fit
- Easy to trim and is compliant with the overlying soft tissues
- Prevents migration of bacteria into wound if exposed.
Edges remain soft and supple to prevent flap complications
- Prevents tissue ingrowth making removal of membrane easier when compared to removal of titanium mesh

Two Different Handling Options

OsseoGuard d-PTFE Titanium-Reinforced is available in two handling options: TR250 or TR150. TR150 membranes are 40% thinner than TR250 membranes, providing clinicians another handling option in titanium-reinforced membranes.



¹ Data on file with manufacturer.


NON-RESORBABLE BARRIER MEMBRANES

OSSEOGUARD PTFE TITANIUM REINFORCED MEMBRANES

Description	Item Number		Units (per box)
	TR250 (250 µm thick)	TR150 (150 µm thick)	
Anterior Extraction 12 mm x 24 mm	TR250AE-1	TR150AE-1	1
	TR250AE-2	TR150AE-2	2
Anterior Extraction 14 mm x 24 mm	TR250AEY-1	TR150AEY-1	1
	TR250AEY-2	TR150AEY-2	2
Large Facial 17 mm x 25 mm	TR250LF-1	TR150LF-1	1
	TR250LF-2	TR150LF-2	2
Posterior Extraction 20 mm x 25 mm	TR250PE-1	TR150PE-1	1
	TR250PE-2	TR150PE-2	2
Posterior 25 mm x 30 mm	TR250P-1	TR150P-1	1
	TR250P-2	TR150P-2	2
Small-T 25 mm x 36 mm	TR250SMT-1	TR150SMT-1	1
	TR250SMT-2	TR150SMT-2	2
Large-T 30 mm x 41 mm	TR250LGT-1	TR150LGT-1	1
	TR250LGT-2	TR150LGT-2	2
Ridge Augmentation X 30 mm x 40 mm	TR250RAX-1	TR150RAX-1	1
	TR250RAX-2	TR150RAX-2	2
Ridge Augmentation K 30 mm x 40 mm	TR250RAK-1	TR150RAK-1	1
	TR250RAK-2	TR150RAK-2	2
Ridge Augmentation K 40 mm x 50 mm	TR250RAKL-1	TR150RAKL-1	1
	TR250RAKL-2	TR150RAKL-2	2
Perio Narrow 13 mm x 19 mm	TR250PN-1	TR150PN-1	1
	TR250PN-2	TR150PN-2	2
Perio Wide 13 mm x 18 mm	TR250PW-1	TR150PW-1	1
	TR250PW-2	TR150PW-2	2
Trans Crestal 24 mm x 38 mm	TR250TCS-1	TR150TCS-1	1
	TR250TCS-2	TR150TCS-2	2
Trans Crestal 38 mm x 38 mm	TR250TCL-1	TR150TCL-1	1
	TR250TCL-2	TR150TCL-2	2
Posterior Ridge 38 mm x 38 mm	TR250PR-1	TR150PR-1	1
	TR250PR-2	TR150PR-2	2



Anterior Extraction
12 x 24 mm




Trans Crestal
38 x 38 mm



Ridge Augmentation K
40 x 50 mm

OsseoGuard® Titanium Mesh

High-Density PTFE Membrane



OsseoGuard Titanium Mesh High-Density PTFE Membranes have similar characteristics as the OsseoGuard Textured, (aside from textured surface). The high-density PTFE membrane is offered in a cost-effective configuration.

Key Attributes¹

- Ultra-thin; 0.2 mm thick
- 0.5 mm pore size
- Highly inert, non-reactive, non-stick nitride coating
- Can be repeatedly sterilized by autoclave

OSSEOGUARD PTFE TITANIUM MESH MEMBRANES


Item Number	Description	Units (per box)
TIM2534-1	Titanium Mesh 25 mm x 34 mm	1
TIM4545-1	Titanium Mesh 45 mm x 45 mm	1

Benefits¹

- 0.5 mm pore size contains graft material while allowing tissue ingrowth
- High coating density with no pores to hold contaminants. Unused portions are not wasted.
- Designed to make primary closure easier to achieve, and to improve tissue release upon removal
- Outstanding wear resistance material will not stain or corrode, and withstands acids, bases, solvents, and high temperatures

OsseoGuard® PTFE

Textured and Non-Textured Membranes



Key Attributes¹

- Non-Resorbable, can be left exposed
- 100% Dense (non-expanded) PTFE
- Soft tissue attaches, but doesn't grow through the membrane

Benefits¹

- Does not resorb prematurely – you dictate healing time
- Impervious to bacteria (pore size less than 0.3 µm)
- Preservation of the soft-tissue architecture and keratinized mucosa
- Exposed membrane allows for non-surgical removal; no anesthesia

OSSEOGUARD PTFE TEXTURED MEMBRANES

Item Number	Description	Units (per box)
TXR1224-1	Textured 12 x 24 mm	1
TXR1224-10		10
TXR2530-1	25 x 30 mm	1
TXR2530-4	25 x 30 mm	4

OSSEOGUARD PTFE NON-TEXTURED MEMBRANES

Item Number	Description	Units (per box)
NTXR1224-10	Non-Textured 12 x 24 mm	10
NTXR2530-4	Non-Textured 25 x 30 mm	4

¹ Data on file with manufacturer.

Collagen Matrices

Plug, Tape, and Patch

Highly porous, absorbable Collagen Wound Dressings to protect, heal, and repair oral wounds.

Clinical Evidence

- Resorption and remodeling profiles of carbonate apatite mimic natural bone mineral¹
- Higher osteoconductive properties and earlier bioresorption, compared to HA samples^{2,3,4}

COLLAGEN MATRICES: PLUG, TAPE AND PATCH

Item Number	Description
0100Z	Zimmer Collagen Tape, 25 x 75 x 1 mm, 10 u/pk
0101Z	Zimmer Collagen Patch, 20 x 40 x 3 mm, 10 u/pk
0102Z	Zimmer Collagen Plug, 10mm x 20 mm, 10 u/pk

Key Attributes

- Made of porcine collagen¹
- Holds up to 30x own weight in fluid²
- No removal needed – resorbs in fewer than 30 days²
- Greater than 90% open pores²
- Protects wound bed – adheres and provides coverage to oral wounds and sores
- Designed to aid healing – porous, absorbable matrix supports delicate new tissue

Clinically successful in procedures for:

- Periodontal surgical wounds¹
- Suture sites¹
- Extraction sites¹
- Surgical wounds¹
- Traumatic wounds¹



Collagen Plug
10 x 20 mm



Collagen Tape
25 x 75 mm, 1 mm thick



Collagen Patch
20 x 40 mm, 3 mm thick

¹ ZimVie Collagen Absorbable Wound Dressings IFU latest revision. ² Data on File with Collagen Matrix Inc.


OsseoGuard®

Non-Resorbable Sutures

Non-Resorbable PTFE Soft Monofilament Suture

Key Attributes

- 100% medical grade PTFE
- Very low package memory
- Biologically inert
- Does not allow bacteria wicking into the surgical site
- Comfortable for patients
- Excellent handling, knots securely
- Keeps the surgical site reliably closed



OSSEOGUARD NON-RESORBABLE SUTURES

Item Number	Description
OS4013PE	USP 4-0, 13 mm, 1/2 circle round body, taper point
OS4013PR	USP 4-0, 13 mm, 3/8 circle precision, reverse cutting
OS3016	USP 3-0, 16 mm, 3/8 circle precision, reverse cutting
OS4016	USP 4-0, 16 mm, 3/8 circle precision, reverse cutting
OS2019	USP 2-0, 19 mm, 3/8 circle precision, reverse cutting
OS3019	USP 3-0, 19 mm, 3/8 circle precision, reverse cutting
OS3016B	USP 3-0, 16 mm, 3/8 circle precision, reverse cutting black
OS3019B	USP 3-0, 19 mm, 3/8 circle precision, reverse cutting black


Safescraper Twist

Cortical Bone Collector

Effectively harvesting autogenous bone which contains viable bone cells which might contribute to the outcome of bone grafting procedures.¹

Key Attributes

- Provides 160° cutting area to effectively harvest up to 5 cc of cortical bone
- Available in curved and straight designs facilitating access to hard-to-reach posterior regions
- Harvested bone is contained in a sterile chamber
- Harvested bone contains viable bone cells and shows high osteogenic potential^{1,2}
- Higher cell viability, cell proliferation, osteogenic potential, and release of growth factors compared to other harvesting methods^{2,3}



SAFESCRAPER TWIST BONE COLLECTOR

Item Number	Description
3598	Disposable Cortical Bone Collector, 3 Units/pk, Straight
3987	Disposable Cortical Bone Collector, 3 Units/pk, Curved


¹Zaffe D. et al. Clin Oral Implants Res (2007) 18:525-533. ²Miron R.J. et al. J Dent Res (2011) 90:1428-33. ³Miron R.J. et al. Clin Implant Dent Relat Res (2013) 15:481-489.

Sinus Crestal Lift Instrument Kit

Crestal sinus lifts without fear of membrane damage.

Key Attributes

- Precisely control drilling depth in increments of 1 mm
- Intuitive layout designed for efficient workflow
- Extended drill sizes with diameter range of 2.4 mm to 4.4 mm



SINUS CRESTAL LIFT KIT

Item Number	Description
SCAKITV2	Sinus Crestal Approach Kit Version 2
SCDRILL24	Sinus Crestal Drill, 2.4 mm
SCDRILL28	Sinus Crestal Drill, 2.8 mm
SCDRILL32	Sinus Crestal Drill, 3.2 mm
SCDRILL36	Sinus Crestal Drill, 3.6 mm
SCDRILL44	Sinus Crestal Drill, 4.4 mm
SCINTDRILL	Sinus Crestal Initial Drill
SCINSERTER	Sinus Crestal Inserter
SCBONESYR	Sinus Crestal Bone Syringe
SCBONECNSDR	Sinus Crestal Bone Condenser
SCDPHTGAUGE	Sinus Crestal Depth Gauge
SCSPRDR27	Sinus Crestal Spreader, 2.7 mm
SCSPRDR31	Sinus Crestal Spreader, 3.1 mm
SCSPRDR39	Sinus Crestal Spreader, 3.9 mm


SINUS CRESTAL LIFT KIT

Item Number	Description
SCSTOPPER02	Sinus Crestal Stopper, 2 mm
SCSTOPPER03	Sinus Crestal Stopper, 3 mm
SCSTOPPER04	Sinus Crestal Stopper, 4 mm
SCSTOPPER05	Sinus Crestal Stopper, 5 mm
SCSTOPPER06	Sinus Crestal Stopper, 6 mm
SCSTOPPER07	Sinus Crestal Stopper, 7 mm
SCSTOPPER08	Sinus Crestal Stopper, 8 mm
SCSTOPPER09	Sinus Crestal Stopper, 9 mm
SCSTOPPER10	Sinus Crestal Stopper, 10 mm
SCSTOPPER11	Sinus Crestal Stopper, 11 mm
SCSTPRHLDR	Sinus Crestal Stopper Holders (10 pack)
SCACASE	Sinus Crestal Kit Case

Sinus Lateral Lift Instrument Kit

A minimally invasive approach to lateral sinus lift.

- Unique angles and tip on the drill allow for faster and more stable drilling compared to previous generation¹
- Three auxiliary slots provided for placing additional drills
- Specially designed drills allow for small lateral window and minimal flap size¹



SINUS LATERAL LIFT KIT

Item Number	Description
SLAKITV3	Sinus Lateral Approach Kit Version 3
SLELV1	Sinus Lateral Elevator 1
SLELV2	Sinus Lateral Elevator 2
SLELV3	Sinus Lateral Elevator 3
SLACASE	Sinus Lateral Kit Case
SLBUR6	Sinus Lateral Number 6 Bur

SINUS LATERAL LIFT KIT

Item Number	Description
SLDR2065	Sinus Lateral Drill, 2mm, 6.5 mmD
SLDR2075	Sinus Lateral Drill, 2mm, 7.5 mmD
SLDR2085	Sinus Lateral Drill, 2mm, 8.5 mmD
SLDR3565	Sinus Lateral Drill, 3.5mm, 6.5 mmD
SLDR3575	Sinus Lateral Drill, 3.5mm, 7.5 mmD
SLDR3585	Sinus Lateral Drill, 3.5 mm, 8.5 mmD

Screw Fixation Instrument Kit

A solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for ridge augmentation procedures.

Key Attributes

- Power grip connection for secure and stable transfer of the screws to the surgical site
- Two color coded systems, Ø 1.5 mm MICRO (blue) and Ø 2.0 mmD MINI (red) screws, for easy and rapid identification of the parts possible and simplifies parts matching
- Modular storage system permits individual configuration
- Autoclavable metal storage tray¹

ASSEMBLED START-UP KIT, ITEM 69.01.10Z

Item Number	Description
69.01.11Z	Tray
75.23.52Z	Screw Driver Handle
75.23.23Z	Screw Driver Insert, Micro, Short
75.23.19Z	Screw Driver Insert, Micro, Long
69.01.09Z	Pilot Drill, Micro, 14mmL
69.01.16Z	Pilot Block Drill, Micro

ADDITIONAL PARTS

Item Number	Description
75.23.21Z	Screw Driver Insert, Mini, Short
75.23.22Z	Screw Driver Insert, Mini, Long
69.01.15Z	Pilot Drill, Mini, Short
69.01.17Z	Pilot Block Drill, Mini

FIXATION SCREWS

Item Number	Description
68.85.83Z	Screws, Micro, 1.5 mmD Self Drilling, 3.5 mmL, 10 pack
68.85.84Z	Screws, Micro, 1.5 mmD Self Drilling, 4 mmL, 10 pack
68.85.85Z	Screws, Micro, 1.5 mmD Self Drilling, 5 mmL, 10 pack
68.85.87Z	Screws, Micro, 1.5 mmD Self Drilling, 7 mmL, 10 pack
68.85.49Z	Screws, Micro, 1.5 mmD Self Tapping, 9 mmL, 10 pack
68.85.51Z	Screws, Micro, 1.5mmD Self Tapping, 11 mmL, 10 pack

Please contact ZimVie for a full list of available replacement parts and optional items.



¹ Screw Fixation System IFU latest revision.

For more information, visit [ZimVie.com](https://www.ZimVie.com)

ZimVie

4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272



Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; and all products are manufactured by one or more of the dental subsidiaries of ZimVie (Biomet 3i, LLC, Zimmer Dental Inc., etc.) and marketed and distributed by ZimVie Dental (formerly Zimmer Biomet Dental) and its authorized marketing partners. Puros Allografts, RegenerOss CC Allografts, and RegenaVate DBM are processed by RTI Surgical, Inc. Some Puros Allografts and CopiOs Pericardium are processed by Tutogen Medical GmbH. Cancell SP is a registered trademark of RTI Surgical, Inc. Tutoplast is a registered trademark of Tutogen Medical GmbH. RegenerOss Cortical and RegenerOss Cancellous Allograft Particulate grafts are distributed by ZimVie and processed by Community Tissue Services. RegenerOss Allograft Putty Plus is manufactured by Interpore Cross International, LLC, a subsidiary of Zimmer Biomet with tissue processed and provided by LifeLink Tissue Bank, an AATB-accredited tissue bank. Endobon is manufactured by Biomet France, Sarl. OsseoGuard PTFE products are manufactured by Osteogenics Biomedical, Inc. IngeniOs products are manufactured by Curasan AG. Safescraper is manufactured by Meta Technologies Srl. Screw Fixation Kits are manufactured by Medicon e.G. Sinus Crestal and Lateral Lift Kits are manufactured by Kaiserprecision Co. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of ZimVie. ZV0040 REV A 05/23 ©2023 ZimVie. All rights reserved.

