



ZIMMER BIOMET

Your progress. Our promise.®




Thoracolumbar Solutions

Trell**ss™ -TS**

Porous Ti Interbody System

Surgical Technique Guide



TrellOss-TS (TLIF Straight) interbodies may be used in various posterior approaches to the anterior spine, and is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and with supplemental fixation.

TABLE OF CONTENTS

Patient Positioning	4
Exposure of Operative Level(s)	4
Decompression	5
Distraction	5
Discectomy and Endplate Preparation	6
Implant Size Selection	7
Implant Preparation and Graft Placement	9
Implant Insertion	10
Implant Removal	11
Implant Sizing	11
Kit Contents	12
Important Information on the TrellOss Porous Ti Interbody System	15

Zimmer Biomet Spine does not practice medicine. Each physician should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training physicians have received.

The following general Surgical Technique Guide is for illustrative purposes only. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as to the best treatment for each patient. Only those individuals with specialized training and experience in spinal surgery should attempt to use the TrellOss Porous Ti Interbody System. Detailed preoperative clinical and diagnostic evaluation followed by carefully executed surgical technique is essential.

Refer to the Instructions for Use (IFU) for a complete list of prescribing information. This technique guide was developed in conjunction with health care professionals.



Figure 1
Patient positioning

PATIENT POSITIONING

- Following adequate general anesthesia, the patient is placed in the prone position on a radiolucent spine table (Figure 1). Particular attention is applied to the positioning of the head and extremities to lessen the risk of ocular and nerve compression.




Figure 2
Confirm correct operative level(s)



Figure 3
Confirm correct operative level(s)

EXPOSURE OF OPERATIVE LEVEL(S)

- Identify the affected level(s) using fluoroscopic imaging and palpation of the targeted anatomy (Figure 2). Access the operative site using preferred instruments. Tissues should be retracted enough to allow for exposure and visualization of the targeted disc space. Insert a marker into the disc(s) to confirm the correct operative level(s) using a lateral radiograph (Figure 3).

Note: TrellOss-TS interbodies are indicated for use at up to two contiguous levels in the lumbar spine, from L2-S1.




Figure 4

DECOMPRESSION


- Utilizing osteotomes and rongeurs, a small section of the lamina and facet(s) should be removed to create an appropriately sized bony window for access to the targeted disc space (Figure 4).




Figure 5

DISTRACTION (optional)

- Effective distraction aids in removal of the superior articular process, decompression of the neuroforamen, preparation of the disc space and insertion of the implant. This may be accomplished through several techniques: pedicle screw distraction, distraction between bony elements, and/or distraction with paddle distractors (Figure 5).



TLIF or PLIF Approach



PLIF – Straight Approach




Figure 7

Figure 6

DISCECTOMY AND ENDPLATE PREPARATION

- Access to the disc space is achieved through an annulotomy made lateral to the posterior longitudinal ligament.
- Using a scalpel, vertical cuts should be made parallel to the dura and laterally in the foramen between the superior and inferior endplate. Additional cuts extend horizontally along the endplates, connecting the vertical cuts.
- Perform a complete discectomy using preferred surgical instruments. Pituitaries, cup curettes (Figure 6), rongeurs and interspace shavers (Figure 7) may be used to remove the disc material. The discectomy is complete once superficial layers of the entire cartilaginous endplates are removed and bleeding bone is exposed. If there is significant disc space collapse, a complete discectomy may not be possible until disc space distraction is accomplished.

- Appropriate endplate preparation will optimize surface contact with the selected TrellOss-TS interbody.

Note: Prior to placement of the implant, autograft or allograft may be placed in the anterior and lateral aspects of the intervertebral disc space.

Warning: Excessive removal of subchondral bone may weaken the vertebral endplate. If the entire endplate is removed, subsidence and a loss of segmental stability may result.




Figure 8



Figure 9

IMPLANT SIZE SELECTION

- Selection of the trial depends on the height, width, and depth of the intervertebral space. Based on the pre-operative imaging and surgical technique, connect an appropriately sized trial to the quick release T-handle and insert it into the annulotomy window (Figure 8).
- Each trial is labeled to differentiate height and should be used incrementally to determine the appropriate dimensions of the required TrellOss-TS implant (Figure 9).




Figure 10




Figure 11

IMPLANT SIZE SELECTION *(continued)*

- Insert desired trial into the intervertebral disc space using gentle impaction of a mallet. Fluoroscopy can assist in confirming the fit and geometry of the trial. If the trial appears too small or too tight, try the next larger or smaller size until the most secure fit is achieved.
- The slap hammer can be used to facilitate the removal of the trial from the intervertebral disc space. The slap hammer includes laser etched steps (attach, drop and rotate) for attachment of the trial (Figure 10). To use, apply an upward force to the slap hammer. Repeat until trial is removed from the intervertebral disc space.

Notes:

Adequate preparation of the endplates is critical in facilitating vascular supply to promote fusion.

Trial sizes (d x w x h) are a line to line match to the corresponding implant. There is no need to under-size or over-size the implant.

All Implant heights are measured from the tallest point on the Implant (Figure 11).

All implants have superior/inferior teeth to help resist implant migration and expulsion while providing a high degree of initial stability.




Figure 12

IMPLANT PREPARATION AND GRAFT PLACEMENT

- Open the sterile packaging of the TrellOss-TS Implant size (height and footprint) that was determined with the trial.
- Prior to insertion, pack the center cavity of the implant with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft. In addition, autograft or allograft may be placed in the anterior and lateral aspects of the intervertebral disc space.



Figure 13



Figure 14

IMPLANT INSERTION

Note: The TrelOss-TS Porous Ti Interbody System comes complete with two inserters specific to the needs of the individual procedure. A minimally invasive (offset) inserter and straight inserter are provided. Both instruments assemble to the implant and function in the same manner.

- Place the end of the Implant into the circular cup and thread into the implant by turning the knob on the threaded shaft component clockwise (Figure 13). Confirm the implant is securely attached but DO NOT overtighten.
- Insert the implant into the intervertebral disc space. If necessary, controlled and light hammering with a mallet can be used to help advance the implant to the desired position (Figure 14).
- The use of fluoroscopy is recommended during any or all of the implantation steps to ensure proper positioning.
- Turn the knob on the threaded inserter shaft in a counterclockwise direction to release the implant from the Inserter (Figure 13).
- If the Implant requires further adjustment, a tamp may be used to carefully manipulate into desired position.
- Complete the procedure by following the surgical technique for the specific device to be used as supplemental fixation.




Figure 15


IMPLANT REMOVAL

- Attach either the inserter or universal removal instrument in a clockwise rotation to the implant (Figure 15). Be careful to avoid pushing the implant anteriorly. A slap hammer or slotted mallet may be used in conjunction with the inserter for removal of the implant if desired. To use, apply an upward force to the slap hammer. Repeat until implant is removed from the intervertebral disc space.
- If distraction was utilized during implantation, be sure to re-apply distraction to allow easier removal of the implant. Vertebral bone overgrowth or osteophytes may be removed to facilitate retrieval of the implant.

Notes:

An osteotome can be used at the interface between the Implant and endplates to disengage the construct.

Use of distraction is suggested to allow easier access to the implant/endplate interface.



IMPLANT SIZING

6° Straight Implants

LENGTHS	ANTERIOR HEIGHT	22 mm	26 mm	30 mm
		POSTERIOR HEIGHT		
8 mm x 6°	8 mm	7 mm	6 mm	6 mm
9 mm x 6°	9 mm	8 mm	7 mm	7 mm
10 mm x 6°	10 mm	9 mm	8 mm	8 mm
11 mm x 6°	11 mm	10 mm	9 mm	9 mm
12 mm x 6°	12 mm	11 mm	10 mm	10 mm
13 mm x 6°	13 mm	12 mm	11 mm	11 mm
14 mm x 6°	14 mm	13 mm	12 mm	12 mm
15 mm x 6°	15 mm	14 mm	13 mm	13 mm
16 mm x 6°	16 mm	15 mm	14 mm	14 mm

KIT CONTENTS

TrellOss-TS 22 mm 0° Implant Kit Kit Number: PCR200M1101

DESCRIPTION	QTY	PART NUMBER
Straight 22D x 9W x 7H 0°	2	207M2207
Straight 22D x 9W x 8H 0°	4	207M2208
Straight 22D x 9W x 9H 0°	4	207M2209
Straight 22D x 9W x 10H 0°	4	207M2210
Straight 22D x 9W x 11H 0°	4	207M2211
Straight 22D x 9W x 12H 0°	4	207M2212
Straight 22D x 9W x 13H 0°	4	207M2213
Straight 22D x 9W x 14H 0°	2	207M2214
Straight 22D x 9W x 15H 0°	2	207M2215
Straight 22D x 9W x 16H 0°	2	207M2216

TrellOss-TS 26 mm 0° Implant Kit Kit Number: PCR200M2101

DESCRIPTION	QTY	PART NUMBER
Straight 26D x 9W x 7H 0°	2	207M2607
Straight 26D x 9W x 8H 0°	4	207M2608
Straight 26D x 9W x 9H 0°	4	207M2609
Straight 26D x 9W x 10H 0°	4	207M2610
Straight 26D x 9W x 11H 0°	4	207M2611
Straight 26D x 9W x 12H 0°	4	207M2612
Straight 26D x 9W x 13H 0°	4	207M2613
Straight 26D x 9W x 14H 0°	2	207M2614
Straight 26D x 9W x 15H 0°	2	207M2615
Straight 26D x 9W x 16H 0°	2	207M2616

TrellOss-TS 30 mm 0° Implant Kit Kit Number: PCR200M3101

DESCRIPTION	QTY	PART NUMBER
Straight 30D x 10W x 7H 0°	2	207M3007
Straight 30D x 10W x 8H 0°	4	207M3008
Straight 30D x 10W x 9H 0°	4	207M3009
Straight 30D x 10W x 10H 0°	4	207M3010
Straight 30D x 10W x 11H 0°	4	207M3011
Straight 30D x 10W x 12H 0°	4	207M3012
Straight 30D x 10W x 13H 0°	4	207M3013
Straight 30D x 10W x 14H 0°	2	207M3014
Straight 30D x 10W x 15H 0°	2	207M3015
Straight 30D x 10W x 16H 0°	2	207M3016

TrellOss-TS 22 mm 6° Implant Kit Kit Number: PCR200M4111

DESCRIPTION	QTY	PART NUMBER
Straight 22D x 9W x 8H 6°	4	208M2208
Straight 22D x 9W x 9H 6°	4	208M2209
Straight 22D x 9W x 10H 6°	4	208M2210
Straight 22D x 9W x 11H 6°	4	208M2211
Straight 22D x 9W x 12H 6°	4	208M2212
Straight 22D x 9W x 13H 6°	4	208M2213
Straight 22D x 9W x 14H 6°	2	208M2214
Straight 22D x 9W x 15H 6°	2	208M2215
Straight 22D x 9W x 16H 6°	2	208M2216

TrellOss-TS 26 mm 6° Implant Kit Kit Number: PCR200M5111

DESCRIPTION	QTY	PART NUMBER
Straight 26D x 9W x 8H 6°	4	208M2608
Straight 26D x 9W x 9H 6°	4	208M2609
Straight 26D x 9W x 10H 6°	4	208M2610
Straight 26D x 9W x 11H 6°	4	208M2611
Straight 26D x 9W x 12H 6°	4	208M2612
Straight 26D x 9W x 13H 6°	4	208M2613
Straight 26D x 9W x 14H 6°	2	208M2614
Straight 26D x 9W x 15H 6°	2	208M2615
Straight 26D x 9W x 16H 6°	2	208M2616

TrellOss-TS 30 mm 6° Implant Kit Kit Number: PCR200M6111

DESCRIPTION	QTY	PART NUMBER
Straight 30D x 10W x 8H 6°	4	208M3008
Straight 30D x 10W x 9H 6°	4	208M3009
Straight 30D x 10W x 10H 6°	4	208M3010
Straight 30D x 10W x 11H 6°	4	208M3011
Straight 30D x 10W x 12H 6°	4	208M3012
Straight 30D x 10W x 13H 6°	4	208M3013
Straight 30D x 10W x 14H 6°	2	208M3014
Straight 30D x 10W x 15H 6°	2	208M3015
Straight 30D x 10W x 16H 6°	2	208M3016

TrellOss-TS Instrument Kit Kit Number: PCR200M4101

DESCRIPTION	QTY	PART NUMBER
Inline Straight Instrument Outer Shaft	1	230M0001
Inline Straight Instrument Inner Shaft	1	230M0002
MIS Straight Instrument Outer Shaft	1	230M0003
MIS Straight Instrument Inner Shaft	1	230M0004
Release Wheel	1	230M0007
Lumbar Removal Tool	1	230M0008
Slap Hammer	1	230M0009
Straight Imp Tamp 90°	1	230M0010
Straight Imp Tamp Flat	1	230M0011
T-handle Silicone Hudson Conn	2	231M0001
Paddle Shaver 6 mm	1	231M0006
Paddle Shaver 7 mm	1	231M0007
Paddle Shaver 8 mm	1	231M0008
Paddle Shaver 9 mm	1	231M0009
Paddle Shaver 10 mm	1	231M0010
Paddle Shaver 11 mm	1	231M0011
Paddle Shaver 12 mm	1	231M0012
Paddle Shaver 13 mm	1	231M0013
Paddle Shaver 14 mm	1	231M0014
Paddle Shaver 15 mm	1	231M0015
Paddle Shaver 16 mm	1	231M0016
Paddle Distractor 6 mm	1	231M1006
Paddle Distractor 7 mm	1	231M1007
Paddle Distractor 8 mm	1	231M1008
Paddle Distractor 9 mm	1	231M1009
Paddle Distractor 10 mm	1	231M1010
Paddle Distractor 11 mm	1	231M1011
Paddle Distractor 12 mm	1	231M1012
Paddle Distractor 13 mm	1	231M1013
Paddle Distractor 14 mm	1	231M1014
Paddle Distractor 15 mm	1	231M1015
Paddle Distractor 16 mm	1	231M1016

TrellOss-TS 0° Trial Kit Kit Number: PCR200M5101

DESCRIPTION	QTY	PART NUMBER
T-handle Silicone Hudson Conn	2	231M0001
Straight Trial 0° 22 x 9 x 7 mm	1	232M2207
Straight Trial 0° 22 x 9 x 8 mm	1	232M2208
Straight Trial 0° 22 x 9 x 9 mm	1	232M2209
Straight Trial 0° 22 x 9 x 10 mm	1	232M2210
Straight Trial 0° 22 x 9 x 11 mm	1	232M2211
Straight Trial 0° 22 x 9 x 12 mm	1	232M2212
Straight Trial 0° 22 x 9 x 13 mm	1	232M2213
Straight Trial 0° 22 x 9 x 14 mm	1	232M2214
Straight Trial 0° 22 x 9 x 15 mm	1	232M2215
Straight Trial 0° 22 x 9 x 16 mm	1	232M2216
Straight Trial 0° 26 x 9 x 7 mm	1	232M2607
Straight Trial 0° 26 x 9 x 8 mm	1	232M2608
Straight Trial 0° 26 x 9 x 9 mm	1	232M2609
Straight Trial 0° 26 x 9 x 10 mm	1	232M2610
Straight Trial 0° 26 x 9 x 11 mm	1	232M2611
Straight Trial 0° 26 x 9 x 12 mm	1	232M2612
Straight Trial 0° 26 x 9 x 13 mm	1	232M2613
Straight Trial 0° 26 x 9 x 14 mm	1	232M2614
Straight Trial 0° 26 x 9 x 15 mm	1	232M2615
Straight Trial 0° 26 x 9 x 16 mm	1	232M2616
Straight Trial 0° 30 x 10 x 7 mm	1	232M3007
Straight Trial 0° 30 x 10 x 8 mm	1	232M3008
Straight Trial 0° 30 x 10 x 9 mm	1	232M3009
Straight Trial 0° 30 x 10 x 10 mm	1	232M3010
Straight Trial 0° 30 x 10 x 11 mm	1	232M3011
Straight Trial 0° 30 x 10 x 12 mm	1	232M3012
Straight Trial 0° 30 x 10 x 13 mm	1	232M3013
Straight Trial 0° 30 x 10 x 14 mm	1	232M3014
Straight Trial 0° 30 x 10 x 15 mm	1	232M3015
Straight Trial 0° 30 x 10 x 16 mm	1	232M3016

KIT CONTENTS *(continued)*

TrelOss-TS 6° Trial Kit Kit Number: PCR200M6101

DESCRIPTION	QTY	PART NUMBER
T-handle Silicone Hudson Connect	2	231M0001
Straight Trial 6° 22 x 9 x 8 mm	1	233M2208
Straight Trial 6° 22 x 9 x 9 mm	1	233M2209
Straight Trial 6° 22 x 9 x 10 mm	1	233M2210
Straight Trial 6° 22 x 9 x 11 mm	1	233M2211
Straight Trial 6° 22 x 9 x 12 mm	1	233M2212
Straight Trial 6° 22 x 9 x 13 mm	1	233M2213
Straight Trial 6° 22 x 9 x 14 mm	1	233M2214
Straight Trial 6° 22 x 9 x 15 mm	1	233M2215
Straight Trial 6° 22 x 9 x 16 mm	1	233M2216
Straight Trial 6° 26 x 9 x 8 mm	1	233M2608
Straight Trial 6° 26 x 9 x 9 mm	1	233M2609
Straight Trial 6° 26 x 9 x 10 mm	1	233M2610
Straight Trial 6° 26 x 9 x 11 mm	1	233M2611
Straight Trial 6° 26 x 9 x 12 mm	1	233M2612
Straight Trial 6° 26 x 9 x 13 mm	1	233M2613
Straight Trial 6° 26 x 9 x 14 mm	1	233M2614
Straight Trial 6° 26 x 9 x 15 mm	1	233M2615
Straight Trial 6° 26 x 9 x 16 mm	1	233M2616
Straight Trial 6° 30 x 10 x 8 mm	1	233M3008
Straight Trial 6° 30 x 10 x 9 mm	1	233M3009
Straight Trial 6° 30 x 10 x 10 mm	1	233M3010
Straight Trial 6° 30 x 10 x 11 mm	1	233M3011
Straight Trial 6° 30 x 10 x 12 mm	1	233M3012
Straight Trial 6° 30 x 10 x 13 mm	1	233M3013
Straight Trial 6° 30 x 10 x 14 mm	1	233M3014
Straight Trial 6° 30 x 10 x 15 mm	1	233M3015
Straight Trial 6° 30 x 10 x 16 mm	1	233M3016

IMPORTANT INFORMATION ON THE TRELLOSS POROUS TI INTERBODY SYSTEM

Device Description

The TrellOss Porous Ti Interbody System is a collection of additively manufactured spacers for cervical, lumbar/lumbosacral and thoracolumbar implantation. The basic shape of these implants is a structural column to provide surgical stabilization of the spine. Each device comprises an external structural frame having a roughened surface (~7 µm). The intervening geometric lattices have pores 300-700 µm. The inferior/superior aspects of the TrellOss open devices incorporate a large vertical cavity which can be packed with bone graft material. The open and solid devices are available in an assortment of height, length, width and lordotic angulation combinations to accommodate the individual anatomic and clinical circumstances of each patient.

Materials

The TrellOss Porous Ti Interbody System implants are manufactured from Titanium Alloy (Ti6Al4V) as described by ASTM F3001.

Indications for Use

- When used as a cervical intervertebral fusion device, the TrellOss-C Porous Ti Interbody System open devices are indicated for use at up to two contiguous levels in the cervical spine, from C2-T1, in skeletally mature patients who have had six weeks of non-operative treatment for the treatment of degenerative disc disease (DDD) with up to Grade 1 spondylolisthesis. DDD is defined as neck pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. The device is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and with supplemental fixation.
- When used as a lumbar intervertebral fusion device, the TrellOss-TS Porous Ti Interbody System open devices are indicated for use at one or two contiguous levels in the lumbar spine, from L2-S1, in skeletally mature patients who have had six months of non-operative treatment for the treatment of degenerative disc disease (DDD) with up to Grade 1 spondylolisthesis. DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. Additionally, the TrellOss Porous Ti Interbody System lumbar implants can be used as an adjunct to fusion in patients diagnosed with degenerative scoliosis. The device is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and with supplemental fixation.

- When used as a vertebral body replacement device, the TrellOss Porous Ti Interbody System open and solid devices are indicated for use in the thoracolumbar spine (T1-L5) for partial replacement (i.e., partial vertebrectomy) of a diseased vertebral body resected or excised for the treatment of tumors or trauma/fracture in order to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body. The device is intended for use with autograft or allograft and with supplemental internal fixation.

Contraindications

The TrellOss Porous Ti Interbody System contraindications include, but are not limited to:

- The presence of infection, pregnancy, metabolic disorders of calcified tissues, grossly distorted anatomy, inadequate tissue coverage, any demonstrated allergy or foreign body sensitivity to any of the implant materials, drugs/alcohol abuse, mental illness, general neurological conditions, immunosuppressive disorders, obesity, patients who are unwilling to restrict activities or follow medical advice, and any condition where the implants interfere with anatomical structures or precludes the benefit of spinal surgery.
- Biological factors such as smoking, use of nonsteroidal anti-inflammatory agents, the use of anticoagulants, etc. all have a negative effect on bony union. Contraindications may be relative or absolute and must be carefully weighed against the patient's entire evaluation.
- Any condition not described in the Indications for Use.
- Prior fusion at the level(s) to be treated.

Warnings and Precautions

- Mixing of dissimilar metals can accelerate the corrosion process. Stainless steel and titanium implants must NOT be used together in building a construct.
- The TrellOss Porous Ti Interbody System devices should be implanted only by surgeons who are fully experienced in the use of such implants and the required specialized spinal surgery techniques. Prior to use, surgeons should be trained in the surgical procedures recommended for use of these devices.
- The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. Based on the dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of the device.

IMPORTANT INFORMATION ON THE TRELLOSS POROUS TI INTERBODY SYSTEM (*continued*)

- The TrelOss solid devices are not intended for interbody fusion as bone growth through the device has not been demonstrated.
- These devices are provided as single use only implants and are not to be reused or reimplanted regardless of an apparent undamaged condition.
- The TrelOss Porous Ti Interbody System is used to augment the development of a spinal fusion by providing temporary stabilization. This device is not intended to be the sole means of spinal support – supplemental internal fixation must be used. If fusion is delayed or does not occur, material fatigue may cause breakage of the implant. Damage to the implant during surgery (i.e., scratches, notches) and loads from weight bearing and activity will affect the implant's longevity.
- The correct handling of the implant is extremely important. Use care in handling and storage of devices. Store the devices in a clean, dry area away from radiation and extreme temperatures and corrosive environments such as moisture, air, etc.
- Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.
- Components of this system should not be used with components of any other system or manufacturer.
- Potential risks identified with the use of this system, which may require additional surgery, include: device component breakage, loss of fixation/loosening, non-union, vertebral fracture, neurologic, vascular or visceral injury.

POTENTIAL ADVERSE EFFECTS

Potential complications and adverse effects for this system are similar to those of other spinal instrumentation systems and include but are not limited to: pseudarthrosis, insufficient bone stock, painful bursa, pressure necrosis, palpable components, early or late loosening of the components; disassembly, bending or breakage of any or all of the components; foreign body (allergic) reaction to the implants; infections possible requiring removal of devices; loss of neurological function, including paralysis, spinal cord impingement or damage.

Reference

1. Mirkovic SR, Schwartz DG, and Glazier KD. 1995. Anatomic Considerations in Lumbar Posterolateral Percutaneous Procedures. Spine 20 (18): 1965-1971.

Disclaimer: This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Rx Only. Please see the product Instructions for Use for a complete listing of the indications, contraindications, precautions, warnings and adverse effects.



Distributed by:

Zimmer Biomet Spine, Inc.
10225 Westmoor Dr.
Westminster, CO 80021 USA
+1 800.447.3625

Manufactured by:

Nexxt Spine, Inc.
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
+317.436.7801
www.NexxtSpine.com



ZIMMER BIOMET
Your progress. Our promise.®

800.447.3625 / zimmerbiomet.com

©2019 Zimmer Biomet Spine, Inc. All rights reserved.

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet Spine, Inc. or one of its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet Spine. This material is intended for health care professionals, the Zimmer Biomet Spine sales force and authorized distributors. Distribution to any other recipient is prohibited.

2523.1-GLBL-en-REV0619