



 **ZimVie** | **MEDIT**

**Your Path to  
Predictable**



## Scanning Made Easy With a Simple Touch

Get the power to take your  
practice to the next level



# Scanning, Now Easier Than Ever

## Medit Plug and Scan

Medit scanners can be directly connected to the PC using a power delivery cable without the need for power hubs and extra cables. This improves mobility and maintenance.




Minimum System Requirement:  
 USB Type-C/USB Power Delivery  
 5V, 3A Power delivery support

# Digital Dentistry With ZimVie



## Digitize Your Treatment With Medit Link




Scan data created with the Medit scanners can be linked to the workflow of various solution providers via Medit Link. The data can also be used with any software that can open standard 3D formats such as STL, PLY, and OBJ.

### What you can do in Medit Link :

Clear aligners | Crown and bridges | Implant-supported restorations | Full and partial dentures | Orthodontics | Implant surgical guides | Sleep apnea devices | Dental models | 3D printing

# Medit Link Fully Open System

Begin your **clinic's data management** with **Medit Link**.







Create and manage scan data matched with patient information in the electronic chart program and link them with treatment solutions of various workflows.

- Scan with free **intelligent features** and **apps** specialized for your dental procedures.
- Improve consultation experience with your patient to **provide better treatment options**.
- Enhance your communication with labs through a **real time cloud-based system**.





Medit Link's cloud storage encrypts and stores data for reliable protection against virus and ransomware attacks. It also enables data sharing and real-time communication between dentists, patients, and lab technicians anytime, anywhere on a PC, a tablet, or any mobile device with an internet connection.

Experience amazing digital dentistry with Medit Scan for Clinics. Medit's cutting-edge software technology makes scanning more convenient for you.

## For scanning assistance

 <p><b>Smart Scan Guide</b></p> <p>Strengthened user guidance provides help for abnormalities detected during scanning.</p>	 <p><b>Smart Scan Filtering</b></p> <p>Enhanced A.I. filtering performance eases the scanning process by filtering out soft tissue such as lips, cheeks, and tongue.</p>	 <p><b>Smart Color Filtering</b></p> <p>Automatically filters registered colors, such as the color of medical gloves, during scanning.</p>	 <p><b>Smart Stitching</b></p> <p>Scan with your own scanning strategy. Fragmented scan data is aligned with automatically detected common area.</p>
--	---	---	---

## For more precise and accurate scan data

 <p><b>Smart Shade Guide</b></p> <p>Automatically distinguishes and suggests the best matching shade from scan data.</p>	 <p><b>A.I. Scan Body Matching</b></p> <p>Secures accurate scan body data in difficult-to-scan areas, such as the posterior region.</p>	 <p><b>A.I. Abutment Matching</b></p> <p>Secures accurate margin data in difficult-to-scan situations.</p>	 <p><b>HD Mode</b></p> <p>Get a sharper scan image in the margin area with the upgraded HD Mode.</p>
--	---	--	--

# Medit i600 Best-in-Class Scanning Experience

## Light

At 245g, it's lightweight with improved grip for comfort.

## Smooth

Optimal imaging speed delivers a steady scan experience.

## Vivid

Transmits high definition scan data via the new optical system used in the Medit i700.

# The Brain of the i700 with a Splash of Color

## Super Accurate, Super Light, Super Smooth.

The digital workflow begins with accurate data. Medit i600 precisely acquires data for effortless digital dentistry.

**\*High accuracy**  
**Full-arch accuracy:  $10.9 \mu\text{m} \pm 0.98$**

\*Acquires accurate scan data from anterior segment to the occlusal surface

# The Simplest. The Easiest.

Take full control of your scanning with the click of a button. No need for numerous cables and hubs! Simply plug and scan using a single cable for direct connection with a PC.



180° Reversible Tip




Scan Start Button



Detachable Cable

# First Look at the Medit i700

Medit *i700* makes the scanning experience a comfortable one for both the dentist and patient. With both powerful hardware and intelligent software, the full-feature Medit *i700* is the key to unleash your clinic's full potential.



# Medit i700 wireless

## Beyond i700's Proven Tech

Still the same **super fast, light, and accurate i700.**  
But now **wireless.**



### No More Wires

Free from wires, scanning is comfortable at any angle. Scan proximal areas effortlessly without constraints in movement.



### Fast, Accurate, Connected

Scanning speed of up to 70 FPS\*. You can experience the amazingly smooth and fast scanning performance of Medit's best-selling model. Enjoy the fast speed of the i700 wirelessly, flawlessly.



### Supercharged Yet Light

Outstanding weight balance even with a battery attached for added comfort.




### Powerful Hardware Meets Versatile Software

Optimized for convenience. Not only from hardware aspects, but also software. i700 wireless, packed with the latest of MEDIT technology.

\*Frames per second: number of frames of scanning data per second

- Full arch accuracy:  $10.9 \mu\text{m} \pm 0.98$
- Up to 70 FPS



### Small Tip

(Sold separately)

For easy application to the posterior region for patients with small oral cavity

### Reversible Tip

No need to twist your wrist for maxillary scans

### Wider Scan Area

Enlarged scanning window captures large areas at once while maintaining the small size of the scanner head

### Auto Calibration

Automatic tool makes calibration a breeze



Auto Calibration Tool (Sold separately)

### Remote Control


Conveniently check scan data from the scanner itself and prevent cross-contamination

### UV-C LED Disinfection

Lower risk of contamination inside the handpiece

### 60 GHz Wireless Technology

For high-speed wireless connection



### Battery

Scan up to 30 cases with 1 battery\*

### Light and Weight Balanced

For stable, comfortable use


When the scanner's cable disappears, so do restrictions to scan movements. Strengthened impact resistance minimizes damage from accidental drops. Weight balanced design of the i700 wireless provides a stable, comfortable scanning experience.

\*Performance may vary according to usage conditions.




# Long-Lasting Battery Life

1 hour of continuous scanning, **up to 8 hours** in standby mode




Scan all day without worrying about charging.



Maxillary + Mandibular + Occlusion  
Total scan time of **2 mins\***




Scan up to **30 cases** with 1 battery



Up to **90 cases** using 3 included batteries



Intelligent power management function automatically switches the device to sleep mode when not in use, standing by up to 8 hours. It takes just seconds to wake the device for scanning.



**Wired power connection:**  
Once you've used up all 3 included batteries, simply connect the power cable and keep using the device.

\*Performance may vary according to usage conditions.



# Medit i900

## The ultimate solution that does it all

Introducing the Medit i900, a premium scanner that redefines the scanning experience - making it not just a task but an effortless and seamless journey for every user.

- Unibody Seamless Handpiece
- Light and Compact
- 360° Touch Band
- Hygienic Touch Pad
- Third Gen New Optical Engine
- Vibrant LED Indicators

## Discover the innovation that sets it apart



**Swipe to Navigate**  
Between Scanning Phases

**360° Touch Band**  
Easy start, stop, and swipe control from any angle

**Touch Pad**  
Effortless navigation through scan data (zoom, pan, rotate)


**Jointless Design**  
Ergonomic comfort and simpler cleaning

**Expanded Field of View**  
Scan larger areas faster

**Detachable Cable**

**Lightweight**  
Just 165 g

**Vibrant LED Light Indicators**  
Clear status updates



# All-New Scanning Experience

\*released

## Third Gen New Optical Engine

- Improved metal and tissue scanning
- Enhanced scan resolution
- Deeper scanning depth
- Better color reproduction



Delivers accurate scan bodies with ease



Provides precise inlay preparation

## Wider Field of View

- 18 mm x 15 mm (552 x 460 px)
- Get precise scans faster with fewer images





Scan reflective surfaces efficiently



Capture detailed full arch scans

# Unleash your clinic's full potential

Category		MEDIT i600	MEDIT i700	MEDIT i700 wireless	MEDIT i900
Model					
Scanning Technology	Capture Rate	Fast	Faster	Faster	Fastest
	Imaging Technology	3D-in-motion video technology/3D full color streaming capture			
	Anti-Fogging Technology	Adaptive air flow			
Accuracy	Full Arch	10.9 um + 0.98			
Handpiece	Dimensions	248 x 44 x 47.4 mm		313 x 44 x 47.4 mm	223 X 37 X 35 mm
	Total Weight	245 g		328 g	165 g
Tip	Scan Area	15 x 13 mm			18 x 15 mm
	Autoclavable	150 times			
		Autoclave 121° C 30 mins (Gravity Type) Autoclave 134° C 4 mins (Pre-Vacuum Type) Autoclave 135° C 10 min (Gravity Type)			Sterilize for 10 minutes at 135°C (275°F) at gravity type and dry for 30 minutes. Sterilize for 4 minutes at 134°C (273°F) at pre-vacuum type and dry for 20 minutes.
Reversible Tip	Yes	Yes	Yes	No	
Special Mode	Remote Control Mode	No	Yes		
	UV-LED Disinfection	No	Yes		
	Smart Shade Guide	Yes			
	HD Mode	No	Yes		
MAC OS Compatibility		Yes			
Connectivity	Connection Type	Wired USB 3.1 Gen 1		Wireless USB 3.1 Gen 1	Wired USB 3.0
	Plug and Scan	Yes			
	Detachable Cable	Yes			

Please scan the QR code for up-to-date [system requirements](#). >>



## About ZimVie

ZimVie is a global life sciences leader in the dental market that develops, manufactures, and delivers a comprehensive portfolio of products and solutions designed to support dental tooth replacement and restoration procedures. From its headquarters in Palm Beach Gardens, Florida, and additional facilities around the globe, ZimVie works to improve smiles, function, and confidence in daily life by offering comprehensive tooth replacement solutions, including trusted dental implants, biomaterials, and digital workflow solutions.

## About Medit

Medit specializes in 3D measurement and CAD/CAM solutions for dental clinics and labs, including intraoral scanners, based on its own patented state-of-the-art technology. We also develop platform solutions for digital dentistry, supporting collaborative workflows. Our goal is to provide innovative technology and the highest quality products to ensure mutual growth for all partners.

## Medit Community

Medit Community provides continuous feedback, suggestions of new treatment workflows, and active discussions via an online community platform. From hot debates to fun chats, the Medit community is always lively with our passionate customers. Colleagues from all over the world actively support and help one another, and sincerely welcome new colleagues into digital dentistry. Scan the QR code to join!

Medit and its community of colleagues support and wish you the best of luck on your new journey into digital dentistry.

**Medit Community:** Connect with users

**Medit Academy:** Check out the resource library where you can learn on the go

**Help Center:** Reach out to a dedicated tech support team

**ZimVie Education:** Check out the world-class ZimVie Institute Programs



Medit User Group



Medit Academy



Medit Help Center



ZimVie Education



[www.medit.com](http://www.medit.com)

[For more information, visit ZimVie.com](http://ZimVie.com)

English

**ZimVie**

4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272



Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; or and affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Intraoral scanners and accessories are manufactured by Medit Corp. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2307 GLBL REV A 06/24 ©2024 ZimVie. All rights reserved.

