

Your path to predictable

Everything you need to take your practice to the next level

Unparalleled technology is here to help you grow your practice efficiency and patient satisfaction:

- Faster, smoother scanning with superior visual output for enhanced metal and edentulous scanning
- Easier handling and advanced ergonomics for more freedom of movement
- Improved scan accuracy for reliable diagnostics
- System includes optimized laptop, support, and a 3-year warranty

Medit i900 scanner

Medit Link

Secure, cloud-based 3D data sharing platform that facilitates communication and collaboration with other dental partners


Compatible with numerous systems, including **RealGUIDE®***, for enhanced workflows.

*anticipated 2024

Medit Apps

A wide range of features to support restorative outcomes:

- Scan and color filtering
- Margin line creation
- Ortho simulation
- Temporaries and crown fit
- A.I. scan body matching



Expanded Field of View
Scan with improved efficiency and speed

360° Touch Band
Easy start, stop, and swipe control from any angle





Touch Pad
Effortless navigation through scan data (zoom, pan, rotate)

Jointless Design
Ergonomic comfort and simpler cleaning

Vibrant LED Light Indicator
Clear status updates



Unleash your clinic's full potential

Category		MEDIT i600	MEDIT i700	MEDIT i700 wireless	MEDIT i900
Model					
Scanning Technology	Capture Rate	Fast	Faster	Faster	Fastest
	Imaging Technology	3D-in-motion video technology/3D full color streaming capture			
	Anti-Fogging Technology	Adaptive air flow			
Accuracy	Full Arch	10.9 um + 0.98			
Handpiece	Dimensions	248 x 44 x 47.4 mm		313 x 44 x 47.4 mm	223 X 37 X 35 mm
	Total Weight	245 g		328 g	165 g
Tip	Scan Area	15 x 13 mm			18 x 15 mm
	Autoclavable	150 times			
		Autoclave 121° C 30 mins (Gravity Type) Autoclave 134° C 4 mins (Pre-Vacuum Type) Autoclave 135° C 10 min (Gravity Type)			
		Sterilize for 10 minutes at 135°C (275°F) at gravity type and dry for 30 minutes. Sterilize for 4 minutes at 134°C (273°F) at pre-vacuum type and dry for 20 minutes.			
Reversible Tip	Yes	Yes	Yes	No	
Special Mode	Remote Control Mode	No	Yes		
	UV-LED Disinfection	No	Yes		
	Smart Shade Guide	Yes			
	HD Mode	No	Yes		
MAC OS Compatibility		Yes			
Connectivity	Connection Type	Wired USB 3.1 Gen 1		Wireless USB 3.1 Gen 1	Wired USB 3.0
	Plug and Scan	Yes			
	Detachable Cable	Yes			

For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272



ZimVie
Education



Medit
User Group



Medit
Academy



Medit
Help Center

 **ZimVie**
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; or and affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Intraoral scanners and accessories are manufactured by Medit Corp. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2397 GLBL REV A 06/24 ©2024 ZimVie. All rights reserved.

